

---

# 2D Multilayer Optics: Construction and Characteristics, Beam Design for Selected Applications and Performance Evaluation

Licai Jiang

Osmic, Inc., 1900 Taylor Road,  
Auburn Hills, Michigan, 48326, USA

E-mail: [licai.jiang@osmic.com](mailto:licai.jiang@osmic.com)



# 2D Multilayer Optical Systems


- Source in point projection
- Graded d-spacing multilayer coating
- Side-by-side optical scheme
- Reflectors of Elliptical or parabolic

Optimized for

- Flux
- Spectral purity
- Spatial definition

Design determining factors


- Source configuration
- Application nature
- Sample characteristics
- Sample size
- Detector characteristics



US patent: 6,041,099

# I would like to talk about these optics

---



CMF



VariMax™




μCMF


# As examples, for these applications

---

Small Molecule  
Crystallography



Protein Crystallography



Small Angle X-ray Scattering


# Regarding these issues

---

- Evaluation of flux
- Evaluation of spatial definition
  - Parallel beam optic: divergence
  - Focusing optic: convergent angle and deviation from the ideal focusing
- Evaluation of spectral purity
  - Continuous background
  - $K_{\beta}$  reduction
- Performance of the systems

# CMF Optics for SAXS

- 3-pinhole setup to eliminate parasitic scattering
- Sample is generally large
- The more difficult experiment is with samples having long periodicity and weak scattering



# CMF Optics for SAXS

---

- Resolution is defined by


$$Q_{\min} = \frac{4p}{l} \sin q_{\min} \quad \text{or} \quad d_{\max} = \frac{l \cdot L_{\text{sample-detector}}}{\left( d_{\text{beam stop}} + d_{\text{beam spot on detector}} \right)}$$

- For the same  $Q_{\min}$ , longer systems yield higher flux
- For the same length, smaller pinholes lead to higher resolution (lower  $Q_{\min}$ )
- Requires a highly defined beam (useful flux is the key)
  - Source with high brilliance
  - Parallel beam optic or focusing beam with low convergence angle (<0.1 degrees)

# CMF Optics for SAXS

---


- Example of a high performance system
  - Rigaku MicroMax007 (70  $\mu$  and 800 W)
  - 2D multilayer focusing optic
  - 3-meter system, > 1000 angstroms configuration
  - more than 15x the flux of a 18 kW rotating anode + graphite monochromator
  - Better resolution





# $\mu$ CMF Optics for SAXS

---

- Another example
  - Bede MicroSource (20x20  $\mu$  and 30 W)
  - 2D multilayer focusing optic  $\mu$ CMF
  - 3-meter system, 1000 angstrom configuration
  - More than 5x the flux of a 18 kW rotating anode + monochromator
  - Better resolution
  - Lower cost (initial and life long)





# $\mu$ CMF Optics for SAXS



Copolymer, with Multiwire Gas Detector

# $\mu$ CMF Optics for SAXS




Turkey tendon, with Multiwire Gas Detector

# $\mu$ CMF Optics for SAXS

---



- Multilayer Wire Gas Detector
  - High sensitivity at low count
  - Real time convenience
- Image Plate Detector
  - High resolution

PET,  
5 minutes exposure with  
image plate



# $\mu$ CMF Optics for SAXS

---




SiO<sub>2</sub>

5 minutes exposure with image plate

RAG setting: 3 kW, 0.2x0.2 mm<sup>2</sup>

# VariMax™ for Single Crystal Diffraction



- Further increased flux over CMF optics
- Adjustable double aperture to change the convergence of the beam
- Improved resolution

# VariMax for Single Crystal Diffraction

---

- Requirement on the beam is more diversified
- Two VariMax have been designed to cover high flux (VariMax HF:  $< 100 - 300$  angstroms); and high resolution (VariMax HR:  $>300 - 500$  angstroms)



- Flux is not always the most important factor
- Beam divergence and intensity (flux/area) plays important role in terms of data quality

# Data Collected from Lysozyme

---

- Size: 0.075 x 0.075 x 0.10 mm<sup>3</sup>
- Optics: VariMax HF and VariMax HR
- Both at high resolution configuration (small apertures)

---

Optic	Range (angstroms)	Ave. counts	R <sub>merge</sub> (shell)	R <sub>merge</sub> (cumul.)
VariMax HF	21.37-1.51	13935	0.058	0.058
VariMax HR	20.61-1.51	5700	0.039	0.039

---

# VariMax Mo coupled with MM007

---



# VariMax for Small Molecule

- Comparing to UltraX 5.4 kW + graphite monochromator (GM)
  - > 10x flux for a sample of 100 micrometers
  - > 4x flux for a sample of 300 micrometers
  - Better beam definition
  - Suppressed high order harmonics






Image at focus

FWHM=0.19 mm


FW10%=0.35 mm

# VariMax for Small Molecule

---



MM007 + VariMax Mo



UltraX 5.4 + GM

# MolyMax™ for Small Molecule

---


- Comparing to 2 kW fine focus sealed tube + graphite monochromator (GM)
  - > 9x flux for 100 micron samples
  - > 3x flux for 200 micron samples
  - > 1.3x flux for 300 micron samples




# MolyMax™ for Small Molecule

---

- $\text{CuC}_{16}\text{H}_{20}\text{N}_2\text{O}_6$
- 0.07x0.11x0.12



Sealed tube



MolyMax

# MolyMax™ for Small Molecule

Source	ST+GM (60 kV, 40 mA)	MolyMax (50 kV, 0.6 mA)
Detector	RAPID	RAPID
Data collection strategy	$c = 45^\circ, f = 0^\circ,$ $0^\circ \leq w \leq 175^\circ$ $c = 45^\circ, f = 90^\circ,$ $0^\circ \leq w \leq 175^\circ$ 5° /image 600 seconds/image	$c = 45^\circ, f = 0^\circ,$ $0^\circ \leq w \leq 175^\circ$ $c = 45^\circ, f = 90^\circ,$ $0^\circ \leq w \leq 175^\circ$ 5° /image 600 seconds/image
Processing package	HKL200	HKL200
Sample: $C_9H_{13}O_5N_3$	$a = 5.1225(7)$	$a = 5.1212(4)$
Size:	$b = 14.0046(17)$	$b = 14.0009(5)$
0.06 x 0.06 x 0.10 mm <sup>3</sup>	$c = 14.7915(11)$	$c = 14.7885(10)$

# MolyMax™ for Small Molecule

---

Source	ST+GM (60 kV, 40 mA)	MolyMax (50 kV, 0.6 mA)
Resolution range (angstroms)	50.0-0.70 (0.73-0.70)	50.0-0.70 (0.73-0.70)
Number of reflections	13316	<b>22305</b>
Average counts	559.3 (30.6)	<b>4226.8</b> (420.3)
$R_{\text{linear}}$	0.066 (0.591)	0.038 (0.234)
$I/\sigma(I)$	30.9 (0.6)	100.2 (9.8)
Completeness (%)	80.8 (13.8)	98.9 (94.2)

# MolyMax™ for Small Molecule

---

Source	ST+GM (60 kV, 40 mA)	MolyMax (50 kV, 0.6 mA)
Refinement Program	SHELXL	SHELXL
R <sub>1</sub> factor [I>2.0σ(I)]:	0.0603 (1591)	0.0345 (2881)
R <sub>w</sub> factor [all data]	0.2449 (2459)	0.0856 (3214)
goodness of fit	1.088	1.052
Number of observations	2459	3214

---

# Some General Comments

# Optimization and design


---

- It is the process to find a beam solution based on the application characteristics
- The performance of the system (solution) is determined by all the components such as source, optic, slits, sample and detector
- The process of the optimization and design is always a trade-off among three fundamental parameters: flux, spatial definition and spectral purity

# Evaluating Flux and Spatial Definition

---

- Flux has to be evaluated with required spatial definition for a certain size of sample
  - Divergence for a parallel beam
  - Convergence angle and deviation from ideal focusing



- Higher the requirement on the low divergence, the lower the flux for a small sample

# Evaluating and Comparing Flux

---

- Different applications require different spatial definition
  - High resolution thin film application requires extremely low divergence that only single crystal monochromators can offer ( $\sim 10$  arc seconds)
  - Small angle x-ray scattering requires a beam with very low divergence ( $< 0.1$  degrees)
  - Protein crystallography can tolerate a beam with moderate divergence ( $\sim 0.1 - 0.25$  degrees)
- Even for the same application, divergence requirement could be different
  - Depends on unit cell size
  - Depends on mosaicity of the sample

# Evaluating and Comparing Flux

---

- Flux is determined by
  - Source configuration: brighter the source, higher the flux with better spatial definition
  - Optical configuration: design for larger capture angle, good source usage efficiency and spatial definition
  - Multilayer configuration
    - heavier element results in higher integrated reflectivity with higher divergence
    - Lighter element has higher peak reflectivity and better spatial definition


# Evaluating Spectrum

---


- Spatial definition is determined by both source configuration and optical geometry
- Bigger the source, worse the spectral purity
- For a focusing beam, spectrum is improved by using a small pinhole near the focus (typically  $K_{\alpha} > 98\%$ )
- For many high resolution applications,
  - only  $K\alpha_1$  is allowed
  - Multilayer is not able to suppress  $K\alpha_1$  and  $K\alpha_2$
  - In this case, a single crystal monochromator is needed

# Improved Spectrum

---



From total reflection  
optic + Ni filter



From multilayer optic

# Customized Optics

---

- For different applications including thin film study, stress analysis, high pressure crystallography, reflective SAXS.....
- For different energies including Cr, Co, Cu, and Mo radiations
- Beam size is not much larger than
  - ~ 4x4 mm<sup>2</sup> for Cr K $\alpha$
  - ~ 3x3 mm<sup>2</sup> for Co K $\alpha$
  - ~ 2x2 mm<sup>2</sup> for Cu K $\alpha$
  - ~ 1x1 mm<sup>2</sup> for Mo K $\alpha$

# Summary

---

- Man-made multilayer reflectors are excellent for beam conditioning for x-ray diffraction
- Such optics feature
  - high flux
  - controlled divergence
  - monochromatic spectrum
- 2D multilayer optics have been applied for
  - Protein crystallography
  - Small molecule crystallography
  - High pressure
  - Small angle x-ray scattering
  - High resolution diffraction
  - Thin-film reflectometry
  - .....

# Acknowledgements

---

I would like to thank

- Dr. Yoshio Iwasaki of Rigaku for the SAXS data collected with image plate detector
- Dr. Scott Barton of Molecule Metrology for the SAXS data collected with Multiwire Gas Detector
- Drs. Joseph D. Ferrara, Lee M. Daniels, and Ron Benson of Rigaku/MSU for the data of protein diffraction and small molecule diffraction