

XRF and SEM-EDS

- Using the PDF for material identification using elemental data from XRF and SEM-EDS.

XRF and SEM-EDS

What?

The Powder Diffraction File contains data on pure solid state compounds of well defined elemental composition.

XRF and SEM-EDS methods can provide an experimental determination of a specimen's elemental composition.

By matching experimental composition data to the database entries, materials can be identified. The more elements positively identified in the XRF or SEM-EDS experiment, the more narrow the selection of candidate phases will be.

If additional data are used, such as physical properties (color, density) or chemical properties, a unique solution can often be identified.

XRF and SEM-EDS

Why?

Scientists have long recognized that using multiple observations of a specimen increases the probability of a successful identification.

The Powder Diffraction File is designed as a database for material identification. While characteristic diffraction and crystallographic data are a primary tool used in the database, other characteristics of a material are input into the database or calculated to increase the chances of a successful identification.

Elemental composition has been experimentally determined or calculated for all entries in the Powder Diffraction File.

Elemental data, from an XRF or SEM-EDS, are often available in global analysis and materials characterization facilities.

XRF and SEM-EDS

How?


All entries in the PDF have calculated atomic and weight percent compositions. The former was designed for use with EDS data and the latter with XRF analyses.

The PDF has composition searches so that experimental data can be compared to the references in the database. ESD's and elemental ranges can be applied to the search.

XRF and SEM-EDS

How?

The more you know, the more efficient the search.



Material Identification

Case 1 - Meteor

A specimen of a commercial meteorite was examined.
The task was to evaluate the composition and verify the authenticity.

What do you know ?

Elemental composition determined by XRF –

Do a composition search.

The specimen is a meteorite –

Use the metal and alloy subfile and/or the mineral subfile.

The specimen was also be examined visually. It appeared metallic, with minor surface corrosion and metallic gray underneath the surface layer. The specimen was heavy.

Data From the XRF Analysis of the Meteor

Compound Formula	nZ	Conc (wt-%)	Stat. Dev. (wt-%)	Line	Peak Int (kCPS)	Bkg Int (kCPS)	Net Int (kCPS)
Na	11	0.27	0.022	Na KA1-HS-Min	0.207	0.035	0.172
Mg	12	0.071	0.0079	Mg KA1-HS-Min	0.250	0.099	0.151
Al	13	0.1	0.011	Al KA1-HS-Min	0.182	0.020	0.162
Si	14	0.23	0.015	Si KA1-HS-Min	0.451	0.025	0.426
P	15	0.47	0.018	P KA1-HS-Min	1.219	0.045	1.174
S	16	0.07	0.0057	S KA1-HS-Min	0.375	0.054	0.321
Cl	17	0.1	0.0099	Cl KA1-HR-Min	0.281	0.047	0.234
K	19	0.042	0.0042	K KA1-HS-Min	0.305	0.062	0.243
Ca	20	0.033	0.0041	Ca KA1-HS-Min	0.292	0.091	0.201
Cr	24	0.017	0.0022	Cr KA1-HS-Min	0.456	0.201	0.255
Fe	26	88.05	0.1	Fe KA1-HS-Min	1167.631	0.514	1167.117
Co	27	0.568	0.0063	Co KA1-HS-Min	14.569	0.556	14.013
Ni	28	9.98	0.1	Ni KA1-HS-Min	15.133	0.074	15.059
Total		100.00					

Composition Search Meteor

Search

Global Operator Numeric Input Help

Subfiles/Database Filters Periodic Table **Elements** Names References Structures Miscellaneous

Empirical Formula
 Not Contains Elements Contains Phrase

Formula Type (ANX)
 Not Contains Exactly

Composition (1)

Not Weight % Atomic %

El.	Value	ESD
Fe	88.0	1.0
Ni	10.0	1.0

Number of Elements (# El's)
 Not Or

1
2
3
4
5
6
7
8
9
≥10

Search Show Results Undock Page Reset Page Reset All

Use the Elements search page in PDF-4+


Composition (2)

Not Weight % Atomic %

El.	Value	ESD
Fe	88.0	1.0
Ni	10.0	1.0

Input major elements from the XRF analysis. Click on a search using weight percents.

Subfile Search Meteor



The screenshot shows the ICDD DDView+ software interface. The main window is titled "ICDD DDView+ - PDF-4+ 2007 RDB". The "Search" window is open, showing various filters and a tree view of subfiles. The "Subfile/Subclass" list includes categories like "Battery Material", "Cement Material", "Ceramics", "Common Phase", "Education", "Expiclient", "Explosive", "Forensic", "Giant Magneto Resistance", "Inorganic", "Intercalate", "Ionic Conductor", "Metals & Alloys", "Mineral Related", "NBS", and "Organic". The "Mineral Related" category is expanded, showing subfiles: "Mineral", "Gem", "Non-ambient", and "Synthetic". A red box highlights the "Mineral" subfile, and a red arrow points from it to a callout box on the right.

Global Operator Numeric Input Help

Subfiles/Database Filters Periodic Table Elements Names References Structures Miscellaneous

Database

Not ICDD (00) Or Not Primary Status Ambient/Non-ambient

Not ICSD-FIZ (01) Or Not Alternate Ambient

Not Cambridge (02) Or Not Deleted Temperature (Non-ambient)

Not NIST (03) Or Not Pressure & Temperature (Non-ambient)

Not LPF (04) Or Not Pressure & Temperature (Non-ambient)

Quality Mark (QM)

Not Star (S) Or

Not Indexed (I) Or

Not Blank (B) Or

Not Low-Precision (O) Or

Not Calculated (C) Or

Not Prototyping (P) Or

Not Rietveld (R) Or

Not Hypothetical (H) Or

Subfile/Subclass

Not Battery Material And

Not Cement Material Or

Not Ceramics Or

Not Common Phase Or

Not Education Or

Not Expiclient Or

Not Explosive Or

Not Forensic Or

Not Giant Magneto Resistance Or

Not Inorganic Or

Not Intercalate Or

Not Ionic Conductor Or

Not Metals & Alloys Or

Not Mineral Related Or

Not Mineral Or

Not Gem Or

Not Non-ambient Or

Not Synthetic Or

Not NBS Or

Not Organic Or

Search Show Results Undock Page Reset Page Reset All

start Windows Expl... Microsoft PowerP... untitled - Paint ICDD DDView+ - ... Microsoft Excel - ... 10:19 AM

Either Mineral or Metal & Alloys, or both, could be selected through the point and click interface.

Mineral

Results

Composition (2)

Not

Weight % Atomic %

El.	Value	ESD
Fe	88.0	1.0
Ni	10.0	1.0

This composition found 0 results when compared to >270,000 entries.

Search broaden by increasing ESD ranges

Composition (2)

Not

Weight % Atomic %

El.	Value	ESD
Fe	88.0	3.0
Ni	10.0	3.0

4 compositions identified. Each reference can then be examined by clicking on the entry.

Results - {Subfile/Subclass (Metals ...

File Edit Fields Results Indexing Help

Results (4 of 272,232)

Search Preference Set: ICDD Defaults

PDF #	QM	Chemical Formula	Compound Name	Common Name	Mineral Name
04-004-8511	I	Fe _{0.9} Ni _{0.1}	Iron Nickel		
04-004-8513	I	Fe _{0.9} Ni _{0.1}	Iron Nickel		
04-008-8473	I	Fe _{0.8782} Ni _{0.1218}	Iron Nickel	α-Fe _{0.8782} Ni _{0.1218}	kamacite
04-008-8475	I	Fe _{0.8782} Ni _{0.1218}	Iron Nickel	γ-Fe _{0.8782} Ni _{0.1218}	taenite


Results

PDF #	QM	Chemical Formula	Compound Name	Common Name	Mineral Name
04-004-8511	I	Fe _{0.9} Ni _{0.1}	Iron Nickel		
04-004-8513	I	Fe _{0.9} Ni _{0.1}	Iron Nickel		
04-008-8473	I	Fe _{0.8782} Ni _{0.1218}	Iron Nickel	α-Fe _{0.8782} Ni _{0.1218}	kamacite
04-008-8475	I	Fe _{0.8782} Ni _{0.1218}	Iron Nickel	γ-Fe _{0.8782} Ni _{0.1218}	taenite

- 2 Entries have a composition of 90% Fe and 10% Nickel and are man-made alloys.
- 2 Entries have a composition of 88% Fe and 12% Nickel, have an identified mineral names, and were found in a meteorite.

The specimen contains 88.1(1)% Fe and 10.0(1)% Ni with many trace elements.

Results



The reference and editors' comment sections for the two mineral phases, Kamacite and Taenite, mention that they are Fe-Ni alloy polymorphs, and were analyzed from the Carlton meteorite in Hamilton, Texas.

The two minerals have different morphologies and crystal habits.

Verification

Meteor

The commercial vendor of the meteorite claimed that the meteor was from the Tambo Quemada meteorite in Peru. Furthermore, it is classified as an Iron, Medium Octahedrite (IIB) meteor with an 8.7% Nickel content.

(Note: This class of meteorites commonly has taenite and kamacite Fe-Ni Minerals. Furthermore, taenites are often Ni rich, even though the reference from the Carlton meteorite was not.)

Once the references point to taenite and kamacite, each of which has a specific crystal habit, the original specimen was reexamined. Visual examination showed evidence of the habit and color described for these minerals. Two different morphologies are clearly observed.

Conclusion

- The commercial vendor claims were verified. XRF data directed the user to specific minerals contained in the PDF database.
- References in the PDF direct the user to cross confirm visual evidence on color and habit.
- The Ni content in the specimen was slightly higher than that of reference Carlton meteorites, and that claimed by the vendors' bulk analysis of the Tambo Quemada meteorite.
- The Ni variation can be easily explained by differences in concentrations of kamacite and taenite, which are both observed in the specimen.

How Searches Reduce the Candidate Lists

>270,000 Entries in PDF-4 Database

957 Entries with Fe & Ni containing compounds

93 Fe-Ni Alloys

24 Fe-Ni Alloy Minerals

2 Fe-Ni Minerals
with 88(2)% Fe

Alternate Search Strategy

Simplier

>270,000 Entries in PDF-4 Database

24 Minerals containing Ni and Fe

3 Fe-Ni Minerals
with 88(2)% Fe

In this case, the third mineral is Haxonite, an Fe-Ni-Co carbide.

This mineral is found in meteorites and has 0.5 wt% Co, and is often mixed with the other 2 minerals.

Since the specimen contains Co, this phase may be also be present.

Mineral Sample

Case 2

In this experiment, the specimen was a commercial raw material intended for a manufacturing process.

The objective was to verify purity and Composition, claimed by the producer.

Mineral Specimen

Case 2

Compound Formula	nZ	Concentration	Line 1	Concentr. 1	Stat. Dev. 1	Depth 1
Al ₂ O ₃	13	56.81	Al KA1-HS-Min	56.81	0.32	29 um
SiO ₂	14	41.5	Si KA1-HS-Min	41.5	0.28	39 um
P ₂ O ₅	15	0.15	P KA1-HS-Min	0.15	0.018	53 um
K ₂ O	19	0.043	K KA1-HS-Min	0.043	0.0064	0.22 mm
CaO	20	0.042	Ca KA1-HS-Min	0.042	0.0071	0.30 mm
TiO ₂	22	1.13	Ti KA1-HS-Min	1.13	0.019	0.53 mm
Cr ₂ O ₃	24	0.038	Cr KA1-HS-Min	0.038	0.0042	0.90 mm
Fe ₂ O ₃	26	0.635	Fe KA1-HS-Min	0.635	0.0071	1.5 mm
SrO	38	0.029	Sr KA1-HS-Min	0.029	0.0036	17 mm
ZrO ₂	40	0.0405	Zr KA1-HS-Min	0.0405	0.0015	24 mm

XRF analysis of a mineral specimen –
all data expressed as oxides.

Convert to Element Concentration

Experimental Data

56.8% Al₂O₃

41.5% SiO₂

Calculated Data*

30.06% Al

19.39% Si

*Can use MW's provided in the ICDD database.

Use a Composition Search

Composition (2)

Not

Weight % Atomic %

El.	Value	ESD
Al	30.1	3.0
Si	19.4	3.0

Input experimental concentrations

This composition only matches 1 entry, and that entry ($\text{Si}_6\text{Al}_{10}\text{O}_{21}\text{N}_4$) is not a Mineral, but a synthetic ceramic!

No Success

Try a New Search

Broaden composition range –
put ESD's at 5 wt%.

Restrict to minerals (eliminates
synthetic ceramics)

New Search

Case 2

5 Minerals identified as candidates

PDF #	QM	Chemical Formula	Compound Name	Mineral Name
00-007-0330	O	$K - Al_4 (Si Al)_8 O_{20} (OH)_4 \cdot x H_2 O$	Potassium Aluminum Silicate Hydroxide	Illite-Montmorillonite, regular
00-012-0625	B	$Mg_2 (Al , Fe)_6 (Si O_4)_4 O_2 (OH)_2$	Magnesium Aluminum Iron Silicate	Yoderite
00-042-0374	O	$K - Na - Al_2 - Si - O - H_2 O$	Potassium Sodium Aluminum Silicate	Erionite-Na, syn
00-047-0356	B	$Na_{39.8} Al_{70.4} Si_{41.2} O_{207.9}$	Sodium Aluminum Silicate	Chabazite-Na, syn
01-074-1394	I	$Mg_3 Al_6 (Si_8 Al_2)_5 O_{21} (OH)$	Magnesium Aluminum Silicate Hydroxide	Kornerupine

None match manufacturer's label.

Verification

The manufacturer's mineral does have several reference compounds in the database. Composition data for these references are shown below.

Database References

Al33.30 O49.37 Si17.33
Al33.30 O49.37 Si17.33
Al33.30 O49.37 Si17.33
Al33.30 O49.37 Si17.33
Al33.30 O49.37 Si17.33
Al33.30 O49.37 Si17.33

Al / Si = 1.92

XRF Analysis

30.1% Al
19.4% Si

Al / Si = 1.55

The Al and Si concentrations, as well as the Al/Si concentration ratio, indicate that the specimen is not pure. The mineral is often found with SiO₂. The lack of additional elements in the XRF would rule out many other minerals. The Al analysis would indicate ~90% purity.

Conclusion

Case 2

- In this case, the analysis indicates that the specimen is not pure.
- Comparison of the data with standards in the database, suggests impurities, which could be verified with some simple additional testing (i.e. light microscopy examination).

XRF Analysis of the Liberty Bell

Copper	64.95–73.10
Tin	24.00–30.16
Lead	1.30–5.47
Zinc	0.25–1.65
Iron	0.00–0.87
Silver	0.14–0.26
Antimony	0.08–0.18
Arsenic	0.19–0.42
Gold	0.02–0.06
Nickel	0.00–0.28



Range reflects 10 specimens taken from the bell in 1960.
Analysis taken from the Liberty Bell Internet site.

Liberty Bell

Case 3

Composition (2)

Not

Weight % Atomic %

El.	Value	ESD
Cu	69.0	5.0
Sn	27.0	4.0

Input Composition (prior slide)

There are 59 references of the Cu-Zn alloy system – bronze.

PDF #	QM ▲	Weight %	Chemical Formula	Compound Name	SYS	Author	Journal
00-034-1092	I	Cu64.00 Mg6.12 Sn29.88	Cu4 Mg Sn	Copper Magnesium Tin	C	Osamura, K., Murakami.	J. Less-Common Met.
04-001-5604	I	Cu64.00 Mg6.12 Sn29.88	Mg Cu4 Sn	Magnesium Copper Tin	C	Osamura K., Murakami Y.	J. Less-Common Met.
00-006-0650	B	Cu64.00 Mg6.12 Sn29.88	Cu4 Mg Sn	Copper Magnesium Tin	C	Gladyszewsky, E. et al.	Dokl. Akad. Nauk SSSR
01-071-0094	B	Cu69.63 Sn30.37	Cu51.68 Sn12.07	Copper Tin	C	Arnberg, L., Jonsson, A., Westman, S.	Acta Chem. Scand., Ser.

Four alloy matches.

Cu64.00 Mg6.12 Sn29.88
Cu64.00 Mg6.12 Sn29.88
Cu64.00 Mg6.12 Sn29.88
Cu69.63 Sn30.37

Liberty Bell

Case 3

Bronze References

04-007-2188	I	Cu61.63 Sn38.37	Cu3 Sn	Copper Tin
04-007-9969	I	Cu61.63 Sn38.37	Cu6.75 Sn2.25	Copper Tin
01-071-7871	I	Cu62.26 Sn37.74	Cu3.02 Sn0.98	Copper Tin
00-026-0564	B	Cu64.09 Sn35.91	Cu10 Sn3	Copper Tin
01-071-7873	B	Cu64.09 Sn35.91	Cu10 Sn3	Copper Tin
04-007-2013	I	Cu64.09 Sn35.91	Cu10 Sn3	Copper Tin
01-071-0122	B	Cu65.50 Sn34.50	Sn11 Cu39	Copper Tin
04-003-1647	P	Cu65.50 Sn34.50	Cu3.12 Sn0.88	Copper Tin
04-007-1295	S	Cu66.29 Sn33.71	Cu40.4 Sn11	Copper Tin
00-031-0485	C	Cu66.34 Sn33.66	Cu40.5 Sn11	Copper Tin
00-031-0486	C	Cu66.34 Sn33.66	Cu81 Sn22	Copper Tin
01-071-0121	B	Cu66.34 Sn33.66	Sn11 Cu40.5	Copper Tin
00-030-0511	C	Cu66.59 Sn33.41	Cu327.92 Sn88.08	Copper Tin
00-030-0510	C	Cu66.62 Sn33.38	Cu41 Sn11	Copper Tin
01-071-7876	S	Cu66.62 Sn33.38	Cu41 Sn11	Copper Tin
01-071-7877	B	Cu66.62 Sn33.38	Cu41 Sn11	Copper Tin
04-004-2878	S	Cu67.91 Sn32.09	Cu41.5 Sn10.5	Copper Tin
01-071-0094	B	Cu69.63 Sn30.37	Cu51.68 Sn12.07	Copper Tin
00-017-0865	B	Cu74.99 Sn25.01	Cu5.6 Sn	Copper Tin
00-031-0487	B	Cu74.99 Sn25.01	Cu5.6 Sn	Copper Tin
01-071-7875	I	Cu75.21 Sn24.79	(Cu0.85 Sn0.15)	Copper Tin
04-007-3939	B	Cu76.68 Sn23.32	Cu0.86 Sn0.14	Copper Tin
04-001-1512	I	Cu82.81 Sn17.19	Cu0.9 Sn0.1	Copper Tin
04-007-7996	I	Cu84.41 Sn15.59	Cu0.91 Sn0.09	Copper Tin
04-003-2428	I	Cu86.03 Sn13.97	Cu0.92 Sn0.08	Copper Tin

There are several bronze phases, (alpha, beta, gamma) of similar composition within the Sn ranges, found in the Liberty Bell. This shows a few of them.

Closest Match



From Cambridge University

Google Search:

<http://www.msm.cam.ac.uk/phase-trans/2005/bell/bell.html>

Metallurgy of Bronze Bells and Castings

H. K. D. H. Bhadeshia

Bronze used for making bells and gongs is essentially an alloy of copper and tin. Copper, containing about 22-24 wt% of tin, is often known as *bell metal*, because it has a pleasing sound quality when struck.

(Note: This citation was not in reference to the Liberty Bell, but elemental analyses of ancient gongs in Korea).

General Searches

Strategies


The prior examples assume that the user has a **quantitative elemental analysis**, so the preferred search mechanism is the **composition search**.

Qualitative analyses can also be used. In these cases, a general **periodic table search** is effective.



Semi-quantitative results might **use combinations** of the composition search with wide ESD's and a periodic table search and/or an elements search.

Any of the above searches can be combined with diffraction data for dramatically improved results – see the Advanced Identification Tutorial for details. The results of any elemental analysis search can be directly fed into the identification Programs, Sleve or Sleve+, as shown in the tutorial.

Periodic Table Search



Click Search Icon From the Toolbar.
Click Periodic Table.

Search

Global Operator Num

Subfiles/Database Filters **Periodic Table** Elements Names References Structures Miscellaneous

	IA	IIA	IIIB	IVB	VB	VIB	VII B	VIII B	IB	IIB	IIIA	IVA	VA	VIA	VIIA	VIIIA		
Period 1	1 H 1.008															2 He 4.003		
Period 2	3 Li 6.941	4 Be 9.012									5 B 10.811	6 C 12.011	7 N 14.007	8 O 15.999	9 F 18.998	10 Ne 20.180		
Period 3	11 Na 22.990	12 Mg 24.305									13 Al 26.982	14 Si 28.086	15 P 30.974	16 S 32.064	17 Cl 35.453	18 Ar 39.948		
Period 4	19 K 39.098	20 Ca 40.08	21 Sc 44.956	22 Ti 47.88	23 V 50.942	24 Cr 51.996	25 Mn 54.938	26 Fe 55.847	27 Co 58.993	28 Ni 58.69	29 Cu 63.546	30 Zn 65.38	31 Ga 69.72	32 Ge 72.61	33 As 74.922	34 Se 78.96	35 Br 79.904	36 Kr 83.80
Period 5	37 Rb 85.47	38 Sr 87.62	39 Y 88.906	40 Zr 91.22	41 Nb 92.906	42 Mo 95.94	43 Tc 98	44 Ru 101.07	45 Rh 102.91	46 Pd 106.4	47 Ag 107.87	48 Cd 112.41	49 In 114.82	50 Sn 118.69	51 Sb 121.75	52 Te 127.60	53 I 126.90	54 Xe 126.90
Period 6	55 Cs 132.91	56 Ba 137.33		72 Hf 178.49	73 Ta 180.95	74 W 183.85	75 Re 186.21	76 Os 190.2	77 Ir 192.2	78 Pt 195.08	79 Au 196.97	80 Hg 200.59	81 Tl 204.38	82 Pb 207.19	83 Bi 208.98	84 Po 209	85 At 210	86 Rn 222
Period 7	87 Fr 223	88 Ra 226.03																
LN:				57 La 138.91	58 Ce 132.91	59 Pr 140.91	60 Nd 144.24	61 Pm 145	62 Sm 150.36	63 Eu 151.97	64 Gd 157.25	65 Tb 158.93	66 Dy 162.5	67 Ho 164.93	68 Er 167.26	69 Tm 168.93	70 Yb 173.04	71 Lu 174.97
AC:				89 Ac 227.03	90 Th 232.04	91 Pa 231.04	92 U 238.03	93 Np 237.05	94 Pu 244	95 Am 243	96 Cm 247	97 Bk 247	98 Cf 251	99 Es 252	100 Fm 257	101 Md 258	102 No 259	103 Lr 260


Inner Operators: Not And Or

Outer Operators: And Or

Search Show Results Undock Page Reset Page Reset All

Periodic Table

Elements Search



The number of elements selected for the search. For example, in Case 1 (meteor), 2 elements account for 98 wt% of the specimen. If 1 element was selected, only single elements would be considered. If 2 elements were selected, all binary alloys would be searched. If both 1 and 2 were selected, all combinations of elements and binary alloys would be searched.

If you have a candidate formula, as in Case 2, it can be searched here.


Composition Search – Shown in all 3 cases.

Case 1 Revisited

Assume Qualitative XRF Data

Input – Fe and Ni into the Periodic Table Search by point and clicking the elements

Input – All 1 and 2 element combinations (elements and binaries)



Case 1 Revisted

Qualitative XRF Data

QM	Weight %	Mineral Name ▲	SPGR
B	Ni100.00	Awaruite	Fm-3m
B	Fe24.08 Ni75.92	Awaruite	Pm-3m
I	Fe24.08 Ni75.92	Awaruite	Pm-3m
B	Fe48.76 Ni51.24	Kamacite	
O	Fe48.76 Ni51.24	Kamacite	P
S	Fe93.71 Ni6.29	Kamacite	Im-3m
B	Fe48.76 Ni51.24	Taenite	
O	Fe82.63 Ni17.37	Taenite	P4132
B	Fe60.82 Ni39.18	Taenite, syn	Fm-3m
B	Fe100.00	Taenite, syn	Fm-3m
S	Fe48.76 Ni51.24	Tetrataenite	Pm
I	Fe48.76 Ni51.24	Tetrataenite	Pm
O	Fe94.34 Ni5.66	Unnamed mineral (NR)	P
B	Fe94.37 Ni5.63	kamacite	Im-3m
B	Fe94.32 Ni5.68	kamacite	Im-3m
B	Fe94.14 Ni5.86	kamacite	Im-3m
B	Fe94.02 Ni5.98	kamacite	Im-3m
B	Fe93.52 Ni6.48	kamacite	Im-3m
I	Fe87.28 Ni12.72	kamacite	Im-3m
I	Fe83.67 Ni16.33	kamacite	Im-3m
I	Fe87.28 Ni12.72	taenite	Fm-3m
I	Fe83.67 Ni16.33	taenite	Fm-3m
I	Fe70.10 Ni29.90	taenite	Fm-3m
B	Fe48.76 Ni51.24	tetrataenite	P4/mmm

You now have 24 candidate materials in 8 mineral families (mineral name)

However, we also know the specimen is claimed to be from a meteorite!

Search for meteorites!


Case 1 Revisited

Qualitative XRF Data

Two good places to search for “meteorites”:

- The **title** of the references
- The editors’ **comments** (includes specimen details)

Reference Search Page



Global Operator Numeric Input Help

Subfiles/Database Filters Periodic Table Elements Names References Structures Miscellaneous

Long Line
 Not Å ESD: Å L1 L2 L3

Strong Line
 Not Å ESD: Å D1 D2 D3

Density
 Not g/cm³ ESD: g/cm³ Dmeas Dcalc Dstruc

Melting Point
 Not ESD: °C °K °F

Color
 Not And Or
 Black Blue Brown Color Missing Colorless Gray Green Metallic


Organic Functional Group
 Not And Or
 >4_Hetero_Atoms_in_Ring(s)
 >5_Fused_Rings
 >9_Membered_Ring
 1_Hetero_atom_in_ring(s)
 1,2_dIOne_O=C-C=O
 2_Fused_Rings
 2_Hetero_Atoms_in_Ring(s)
 3_Fused_Rings

White Yellow

Smith-Snyder Figure Of Merit (SS/FOM)
 Not ESD:

Database Comments
 Not Contains Words Contains Phrase

Search Show Results Undock Page Reset Page Reset All



Search

Global Operator Numeric Input Help

Subfiles/Database Filters Periodic Table Elements Names References Structures Miscellaneous

Author
 Not Contains Words Contains Phrase

Coden
 Not Contains Exactly

Journal
 Not Contains Words Contains Phrase

Title
 Not Contains Words Contains Phrase

Chemical Abstracts Registry Service (CAS)
 Not Contains Exactly Year
 Not ESD:

Search Show Results Undock Page Reset Page Reset All

Case 1 Revisited

Qualitative XRF Data

- Searching “meteor” in the comment section finds 18 entries of Fe-Ni composition – and four minerals.
- Searching “meteor” in the title finds 13 entries of Fe-Ni composition and three mineral types – including taenite and kamacite.

Conclusions

Since the PDF is a collection of pure single phase materials, identification is enhanced when the specimen is phase pure. This means that this application improves with either XRF or SEM-EDS microanalysis.

The literature citations, reference histories, physical and chemical properties in the database, can all be used to cross reference with experimental elemental data to assist in material identification.



Thank you for viewing our tutorial.

Additional tutorials are available at the ICDD web site
(www.icdd.com).

International Centre for Diffraction Data

12 Campus Boulevard

Newtown Square, PA 19073

Phone: 610.325.9814

Fax: 610.325.9823