

In situ Measurement and Characterization of Crystal Growth by X-ray Diffraction

Crystallization Monitoring – a QbD tool

<u>Detlef Beckers,</u> Klaus Bethke, PANalytical B.V. The Netherlands

This document was presented at PPXRD -Pharmaceutical Powder X-ray Diffraction Symposium

Sponsored by The International Centre for Diffraction Data

This presentation is provided by the International Centre for Diffraction Data in cooperation with the authors and presenters of the PPXRD symposia for the express purpose of educating the scientific community.

All copyrights for the presentation are retained by the original authors.

The ICDD has received permission from the authors to post this material on our website and make the material available for viewing. Usage is restricted for the purposes of education and scientific research.



PPXRD Website – <u>www.icdd.com/ppxrd</u>

ICDD Website - www.icdd.com



Contents of presentation

- Motivation
- Experimental set-up
- *In situ* monitoring of DL-Alanine crystal growth
- Crystal morphology investigation by CT



Crystallization monitoring: a key industry requirement

- Crystallization is always a key aspect of pharmaceutical manufacturing and development:
 - Significant impact on the efficiency and profitability of the overall process
 - Over 90% of all pharmaceutical products contain active ingredients produced in crystalline form
 - Undetected fluctuations in the crystallization process can alter the crystal structure or stability, affecting the safety and the bioavailability of the product.
 - Failure to meet product specifications incurs significant costs



Industry need for on-line monitoring

- Monitoring of crystallization processes both in research and scale-up is essential for a QbD (Quality by design) approach and to develop "PAT" solutions
 - Process Analytical Technology (PAT) is the design and control of manufacturing processes through real-time measurements with the goal of ensuring final product quality
- Paramaters to control: crystal size distributions, polymorph formation and intermediates, impurity-crystal interactions, morphology (shape)
- Current monitor techniques:
 - NIR, Raman: polymorphs
 - FTIR: solution concentration



Experimental – set-up





Applications for *in situ* studies

- Crystallization from supersaturated solutions: e.g. investigation of intermediates (polymorphs) and hemi-hydrates during the crystallization process
- Solvent / anti-solvent reactions (e.g. in cleaning processes , re-crystallization)
- Parameter optimization in the crystallization process (pH, (anti-)solvent concentration, etc...)
- Small angle X-ray scattering (SAXS) studies of early crystallization stages or nano particles
- Scale-up investigation



Sensitivity test - Lactose in ethanol





Experimental: DL-Alanine crystallization

 24.5 g DL-Alanine (C₃H₇NO₂) dissolved in 100ml water at 56°C and cooled down



















12













15



DL-Alanine crystallization at pH = 6.1 - 2D data



16

X-ray tomography - principle

DL-Alanine: different crystal shapes ("morphology")

Conclusions 1 – DL-Alanine crystrallization

- Different crystallization conditions show distinct differences in crystallization initiation and crystal morphology.
- At pH = 6.1 and pH = 9.5, first crystallization is most pronounced by the (002) and (311) reflections. Only much later followed by fast growth of the (210) reflection.

Conclusions 2 - DL-Alanine crystrallization

- Crystallization at pH = 6.1 solution starts after 13 h, whereas at pH = 9.5 first crystals are formed much faster (4 h).
 At pH = 3.5 crystallization is very slow (start after 33 h).
- pH = 6.1; pH = 9.5: peaks become much sharper in the course of the crystallization, indicating increasingly larger crystals.
- In the slow crystallizing condition at pH = 3.5, the crystal morphology is less pronounced - broad peaks point towards smaller crystallites.

Outlook - Ludox® TM-50 silica nanoparticles

