

Study of Biological Macromolecules with a Laboratory XRD system

<u>D. Beckers</u>, T. Degen, G. Nenert, J. Bolze, PANalytical B.V., Almelo The Netherlands S. Logotheti, S. Saslis, I. Margiolaki, University of Patras, Greece V. Kogan, Dannalab, Enschede, The Netherlands

This document was presented at PPXRD -Pharmaceutical Powder X-ray Diffraction Symposium

Sponsored by The International Centre for Diffraction Data

This presentation is provided by the International Centre for Diffraction Data in cooperation with the authors and presenters of the PPXRD symposia for the express purpose of educating the scientific community.

All copyrights for the presentation are retained by the original authors.

The ICDD has received permission from the authors to post this material on our website and make the material available for viewing. Usage is restricted for the purposes of education and scientific research.



PPXRD Website – <u>www.icdd.com/ppxrd</u>

ICDD Website - www.icdd.com

Protein analysis



- Challenges
- Protein poly-crystallography
- In situ humidity studies
- SAXS on bio-macromolecules

The challenge



Proteins are challenging samples:

- Weak scatter
 - \rightarrow high intensity required (and low background)
 - \rightarrow linear detector / area detector (with high resolution)
- Large molecules / cells
 - → good low angle performance (peak position and asymmetry/ resolution)
 - \rightarrow high angular resolution
- Often not stable under radiation



HEWL – Sample in capillary



- HEW Lysozyme
 - Can easily be extracted
 - Well-known procedure for crystallization
 - Crystallization after approx. 46 hours
- Pipetting into 0.5 mm capillary (manually compacted)



Tetragonal

Acknowledgement:

sample supplied by B. Prugovečki, Laboratory of General and Inorganic Chemistry, Faculty of Science, University of Zagreb, Croatia



Intensities drop during ~83 hours radiation (5x17 hours)



HEWL Lysozyme – capillary data





Tetragonal - Orthorhombic



© Copyright (2015), PANalytical B.V.

HighScore Plus

PANalytical

get insight

Bovine Insulin on HTS stage

• Crystallized at ph 5.6

2000

Intensity (counts)

Placed into a HTS-plate





© Copyright (2015), PANalytical B.V.

Acknowledgement: sample supplied by I. Margiolaki, ESRF, Grenoble, France

Indexing





Protein polycrystallography



Insulin - Pawley fit, 4563 refined parameters, 5457 data points



Protein in situ humidity study



- Protein molecules in microcrystalline precipitates are surrounded by solvent and their packing arrangement is retained by limited intermolecular contacts.
- A change in the crystal environment first affects the bulk solvent that fills the intermolecular space, with resulting changes in the crystal structure.

Protein in situ humidity study



Sample preparation:

- Crystallization (HEWL with 50 mg/ml, pH 4.4)
- Concentration in centrifuge
- Pipetting concentrate onto transmission holders (~100µl precipitate per sample) – sample height varies during humidity cycle
- Multiple samples to reduce radiation damage
- Humidity variation: $95\% \rightarrow 51\%$





Protein in situ humidity study



- 3 hours / humidity step
- 30 min stabilization / 2.5 h data collection (3 scans)
- Here: results of 2nd humidity cycle





9.8

94 92 90 Ln(Counts) Decreasing rel. humidity 88 10.33 10.302 86 10.274 10.246 84 10.218 10.19 82 10.162 10.135 80 10.107 10.079 78 10.051 10.023 Humidity [%] 76 9.995 9.967 74 9.939 72 9.911 9.883 70 9.856 9.828 68 9.772 66 9.744 9.716 64 9.688 9.66 62 9.632 9.605 60 9.577 9.549 58 9.521 9.493 56 9.465 54 52 2.8 3.4 3.6 3.8 4.2 5.2 5.4 5.8 3.2 4.4 4.6 4.8 5 5.6 3 4 Position [°20] (Cu K-a12)

Surface view showing phase transitions



- Tetragonal HEWL phase at 95% Humidity
- Pawley fit





Pawley fit of Tetragonal phase from 95% to 79% rH





Peaks of an un-indexed HEWL polymorph after Pawley fit of the Tetragonal HEWL cell at 75% RH





Peaks could be explained by a second tetragonal phases



[©] Copyright (2015), PANalytical B.V.



Pawley fit with two Tetragonal phases below 79% rH



BioSAXS



The biomacromolecules are studied in their native environment, i.e. in dilute (approx. 1 wt.%) solution. Crystallization is not required. Available sample quantities are often very small.

🗸 Size

- Shape (3D envelope)
- Compactness and aggregation behavior
- Folding / unfolding
- 2nd virial coefficient (molecule interaction)

These properties can be studied as a function of temperature, pH, buffer composition, etc.



Glucose isomerase



- Enzyme (protein) produced by microorganisms (bacteria)
- Catalyzes the conversion of glucose to fructose
- Massive industrial use for the production of high-fructose syrups
- Popular protein to demonstrate SAXS measurement and analysis capabilities





Guinier plot / Pair distance distribution function p(r)



	<i>R_g</i> from Guinier plot	R _g from p(r)	D _{max} from p(r)
Lab system	3.31 nm	3.32 nm	9.7 nm
Synchrotron (lit)	3.25 nm	n.a.	9.7 nm

Kratky plot





3D-shape reconstruction





Software: Dammin (EMBL)¹

Model fits the experimental data very well

¹ D. I. Svergun (1999). Restoring low resolution structure of biological macromolecules from solution scattering using simulated annealing. *Biophys J.* 2879-2886.



Comparison with results from single crystal diffraction



Simulation of SAXS data from known single crystal structures (using CRYSOL¹ from the EMBL): **Experimental SAXS** data is in good agreement with the Tetramer structure. (R_g values from simulations)

¹Svergun D.I., Barberato C. and Koch M.H.J. (1995) CRYSOL - a Program to Evaluate X-ray Solution Scattering of Biological Macromolecules from Atomic Coordinates *J. Appl. Cryst.*, **28, 768-773.**



Multi-purpose X-ray laboratory diffraction systems offer a range of possibilities to study biological macro molecules:

- Protein polycrystallography (indexing, space group determination)
- Polymorph identification / screening
- Analysis of structural changes of protein precipitates upon dehydration and hydration
- Shape analysis, aggregation and folding of diluted proteins in solution (SAXS).