

#### T. Hartmann

STOE & Cie GmbH Darmstadt, Germany Email: hartmann@stoe.com

#### **PDF** calculation

# OPTIMISING PDF DATA QUALITY USING A LABORATORY POWDER DIFFRACTOMETER

An impressive comparison of G(r) calculated with PDFgetX3 [1] from data of naphthalene taken at room temperature with a STOE STADI P powder diffractometer in Transmission mode, equipped with a Ag-tube, a Ge(111)-monochromator for pure Ag  $K_{\alpha 1}$ -radiation as well as the Dectris MYTHEN 1K detector with 1 mm chip thickness and from the synchrotron beam line X17A, NSLS Brookhaven with a wavelength of 0.1839 Å resulted into amazingly similar peak widths for both experiments.

Therefore the same powder diffractometer has been chosen to investigate the influence of the wavelength, sample volume, temperature and detector technique to optimise the data quality for laboratory PDF data acquisition.

In the recent past, the calculation of the pair distribution

#### **Temperature dependent measurements**



function (PDF) became more and more common as a new method to obtain structural information from complex materials beyond the Bragg equation and crystallographic PXRD analysis [2].

The atomic PDF G(r) gives the number of atoms in a spherical shell of a normalised thickness at a distance r from an atom,

STADI P with Dectris MYTHEN 1K detector

 $\mathbf{G}(\mathbf{r}) = 4\pi \mathbf{r}[\mathbf{p}(\mathbf{r}) - \mathbf{p}_0]$ 

it is the sine Fourier transform of the total scattering structure function S(Q)

## $G(r) = \frac{2}{\pi} \int_{-\infty}^{\infty} Q[S(Q) - 1]\sin(Qr)dQ$

where Q is the magnitude of the scattering vector and S(Q) the corrected and normalised powder diffraction pattern of the material. S(Q) has to be measured to as high Q-values as possible to achieve the best resolution in real-space.

## The influence of the wavelength

With

 $|Q| = \frac{4\pi \sin \theta}{2}$ 

it is evident that synchrotron radiation with shortest wavelengths will be preferred to attain the best PDF data. Nevertheless, standard laboratory powder diffractometers can yield sufficient data quality in a reasonable measuring time when Mo  $K_{\alpha 1}$ - (0.7093Å) or Ag  $K_{\alpha 1}$ -radiation (0.5594Å) has been used.

Sealed in a 1.0mm capillary, the sample has been measured with 3 different wavelengths. PDF calculations yield a  $Q_{(obs max)}$  for Cu- of 7.0 Å<sup>-1</sup>, 11.4 Å<sup>-1</sup> for Mo-, 13.1 Å<sup>-1</sup> for Ag-K<sub> $\alpha 1$ </sub>-

To observe the temperature dependence of the data quality, the STOE STADI P has been equipped with an additional Oxford Cryosystems Cobra and a STOE furnace

### **Naphthalene at low-temperature**



Room-temperature (red) and -100°C (blue) measurement of naphthalene measured with Ag- $K_{\alpha 1}$ -radiation

hkl	2 <del>0</del> (RT)	2 <del>0</del> (-100°C)	FWHM (RT)	FWHM (- 100°C)
001	4.374	4.444	0.039	0.040
110	6.869	6.894	0.032	0.031
-111	7.053	7.155	0.042	0.041
-201	7.771	7.887	0.053	0.051

The low-temperature experiment of naphthalene has been carried out at -100°C. The comparison with the measurement at room temperature shows the expected shift to higher  $2\theta$ -values, yet the influence on the FWHM is nearly neglectable.



G(r)-plot of naphthalene at RT (red) and -100°C (blue)

## **Ammonium Nitrate at high temperatures**

radiation and 19.5 Å<sup>-1</sup> for the synchrotron data ( $\lambda = 0.1839$  Å).



G(r) for naphthalen measured with Cu- (red), Mo- (green) and Ag-radiation (blue)

#### Synchrotron versus sealed tube



The direct comparison of the PDF curves of the synchrotron- (yellow) and the Agexperiment (blue), shows that the resolution of the Ag-data is amazingly similar!

Taking into account that  $\lambda_{(synchrotron)}$  has been approximately 50% of Ag-K<sub> $\alpha 1$ </sub>, the measuring time (Ag-experiment 18h, synchrotron 10min) is more than reasonable for this laboratory setup making the STOE STADI P with Ag-tube and Dectris MYTHEN 1K an impressive alternative to the expensive synchrotron experiments.

With its 4 phase transitions at moderate temperatures NH<sub>4</sub>NO<sub>3</sub> seemed to be the right



Guinier plot of a rapid T-resolved measurement of NH<sub>4</sub>NO<sub>3</sub>

The comparison of the PDF curves of all 4 phases shows a similar slope up to about 4 Å. Only the 6-fold coordination of N2 yields a remarkable high intensity at r = 2.2 Å in the PDF curve of phase I. For r values greater than 4 Å, the different structural details of the polymorphs are represented in the PDF curves.

Polymorph	N1 – H [Å]	N2 – O [Å]	N1 – O [Å]	N1 – N2 [Å]
Ι		2.2	3.1	3.8





Comparison of the laboratory setup (blue) and the synchrotron experiment (yellow)

#### The right detector technique



Background and powder pattern of PY213 measured with Cu- and lin. PSD, Mo and IP PSD as well as Agradiation and MYTHEN 1K

Pigment Yellow 213 [3] has been chosen as the test sample to compare different detector types, a STOE IP-PSD imaging plate (below), STOE's linear PSD (middle) and the Dectris MYTHEN 1K (top).

Regardless the wavelengths they have been measured with, only a short view the background of the powder on diffraction patterns of pigment yellow 213 (blue) and the empty capillary (red) is enough to conclude that the MYTHEN 1K should be the detector of choice!

II	1.0	$1.2 - 1.25^*$	2.9-3.1*	3.7-3.8*
III		1.3	2.8-2.9	2.6
IV	1.0	1.3	3.0	3.5

G(r)-plot of the  $NH_4NO_3$  polymorphs I - IV

Selected distances for all crystalline  $NH_4NO_3$ -Phases (ICSD data) \*: disorder

#### ACKNOWLEDGEMENTS

We thank Dr. Lothar Fink, Institute of Inorganic and Analytical Chemistry, Goethe-University Frankfurt, Germany and Dr. Nadine Rademacher, Institute of Geosciences, Goethe-University Frankfurt, Germany for the PDF calculations, the group of Prof. Dr. Simon Billinge for the measurements at NSLS, Brookhaven, and Dectris LTD in Baden, Switzerland for the good cooperation.

#### REFERENCES

- [1] Juhas, P., Davis, T., Farrow, C.L. and Billinge, S.J. L., J. Appl. Cryst. 2013, 46, 560-566.
- [2] Dinnebier, R. E. and Billinge, S. J. L., "Powder Diffraction: Theory and Practice", 2008, The Royal Society of Chemistry, Cambridge, UK.
- [3] Schmidt, M.U. et al., *Acta Cryst.* **2009**. *B65*, 189-199.

WWW.STOE.COM - YOUR PARTNER IN X-RAY DIFFRACTION

# This document was presented at PPXRD -Pharmaceutical Powder X-ray Diffraction Symposium

Sponsored by The International Centre for Diffraction Data

This presentation is provided by the International Centre for Diffraction Data in cooperation with the authors and presenters of the PPXRD symposia for the express purpose of educating the scientific community.

All copyrights for the presentation are retained by the original authors.

The ICDD has received permission from the authors to post this material on our website and make the material available for viewing. Usage is restricted for the purposes of education and scientific research.



PPXRD Website – <u>www.icdd.com/ppxrd</u>

ICDD Website - www.icdd.com