ADVANCED MINERAL ANALYSES BY X-RAY POWDER DIFFRACTION

T. G. Fawcett, S. N. Kabekkodu, J. Blanton, C. E. Crowder International Centre for Diffraction Data, Newtown Square, PA, 19373

Acknowledgement

This tutorial was made possible by the many contributions of ICDD members working in the fields of mineralogy and geology.

The data shown in this presentation were contributed by ICDD members and/or clinic instructors

Whats new in X-ray diffraction

- Miniaturized sources
- Energy selective detectors
- Low noise, high dynamic range detectors
- Sophisticated self-aligning optics
- More 2D detectors (texture, particle distributions)
- Improved efficiency multilayer optics
- Growth of worldwide synchrotron facilities and their user communities

Enabling New Instruments for exploration and analysis





XRD and XRF on Mars







Powerful Personal Computers – *Better Analysis Programs*

Rietveld Refinement is commonplace
Whole pattern fitting methods (LeRail)

- Whole pattern fitting methods (LeBail, FullPAT, FULLPROF, GSAS, PONCKS)
- Pair Distribution Function Analyses
- Modulated structure analyses
- Complex multivariant analyses

ICDD Global Users Survey - 2006, 2009 and 2012



Better Data and analysis tools



- Suanite.int — Mg2 B2 O5 - 04-009-3360 (Calc, Intensity: 80.1%) — Mg Fe2 O4 - 04-010-6157 (Calc, Intensity: 22.0%) — C - 04-016-4291 (Calc, Intensity: 55.0%) - Mg B O2 (O H) - 00-019-0755 (Exp-based, Intensity: 4.0%) — Fe - 04-003-5611 (Calc, Intensity: 24.8%)

Instrumental Peak Resolution LaB6 – Better optics and detectors

INCIDENT MONOCHROMATORS

0.30 0.28 LaB6 0.26 GUINIER 0.24 1988 Round Robin 0.22 Powder Diffraction, 3, pp 209–218 0.20 0.1 н 0.1 1 ο, 0.1 0.0 0.0 120 110 DEGREES TWO THETA 0.04 0.026 - Lab 2005 0.02 0.003 – Synchrotron 0.00

Experimental Data – LaB6, 2009 Brian Toby, Argonne National Light Source



9

Laboratory Data



10

Star Quality (Highest quality) Mineral Data in the PDF by year

XRD is a "fingerprint" technique the specimen is always compared to a reference. The highest quality reference data have been produced with modern methods and instrumentation.



41 % of the highest quality data are new since the year 2000

Whats new in XRD analyses of minerals

- Higher quality reference data
- More data
- Better experimental data (optics, sources)
- High efficiency, energy selective, position sensitive detectors
- Improved methods of analysis

Database – Growth with growth in technology

Growth of quality data references based on better experimental data

All entries expressed as digital patterns to enable whole pattern fitting methods

All calculated digital patterns can be adjusted for instrument function and crystallite size including X-ray, neutron, electron and synchrotron diffraction – This enables the study of nanomaterials.

Addition of non-crystalline reference materials experimental patterns - clays

Addition of data mining capability and software for improved phase identification and quantitative analysis by multivariate analysis

Data Mining

Composition

	8A.	HA	-	IVB	VB	VIB	VIIB		VIIIB		18	IIB	IIIA	IVA	VA	VIA	VIIA	VIIIA
-	H		Booler	Boolean Yes/No/Maybe											2 He			
	L	Be		Inter Operators Duter Operators							* 8	C	N	0	*F	10 Ne 20 160		
ariad 3	Na	Mg	E	And Date Own And Or OAnd				nd O	or	13 Al 211 Mil	14 Si 31.000	P	10 S	CI In sta	Ar			
*	19 K	20 Ca 40.075	21 Sc 44.300	12 11 12 M	23 V 10.341	24 Cr 51.880	25 Mn 14.335	20 Fe 85.540	27 Co NJ. HID	28 Ni 14.893	29 Cu 61.748	20 Žn 61.408	31 Ga 55.723	92 Ge 12.64	10 As 74.522	Se TLM	35 Br 79.808	M Kr al.788
2	27 Rb (1.45)	Sr Frid	28 ¥ 86.508	*0 Zr 91.224	Nb SU ND	Mo No	Tc (M)	Ru	Rh	Pd	Ag	Cd	In	50 Sn 118.71	51 Sb 121.78	Te Te	1 1 1 1 1 1 1	54 Xe 191.2NE
	68. CS 132 908	Ba		72 Hf 175.40	73 Ta 101.340	W	Re	70 Os 18121	77 Ir 102,217	Pt Het STS	Au	10 Hg	n 1	Pb	83 Bi 226.98	Po	At [24]	80 Rn (222)
T.	Fr pr	Ra																
.N			E7 La 135.306	Ce	Pr 142.005	Nd 144.242	Pm [141]	82 Sm 182.30	Eu	Gd	TD TD	00 100 100	Ho	68 Er 107.200	Tm 100.304	70 Yb 11104	Lu	
ĸ			Ac III	Th	Pa 211.000	N2 U 200.029	Np (III)	Pu pu	M Am [24]	64 Cm [247]	Bk perj	Cf (251)	88 Es (212)	100 Fm (257)	101 Md (254)	102 No (200)	103 Lr (200)	
				1	Search	14 5	bove Re	auts [/ Une	dock Pag		Reset F	age	Rese	LAL			

Physical Properties

ubfiles/Database Filt	ters Periodic Table Elements Namer	References Structures Miscellaneous	
Long Line		Strong Line	
-	100000		01 -
INCE	A ESO:	A 2 -	A ESD: A 02
Density		1/2-corundum (1/2c)	
		Omean .	
Not	g/cm1 ESD: g/o	m ¹ Deale	ESO:
Malting Print		Briathy	
Not	ESD: O PC	CAK OFF	(50)
Color		Orbanic Functional Group	
Not Black	* O And	Next >4_Hetero_stoms_st	n_ring(s)
Blue		>5_fused_rings	
Brown		>9_membered_ring	dearle)
Color Mesar	*	1_Hetero_etom_in_r	C=0 *
3.00.00			10 ×
Smith-Snuder Finure	e Of Merit (SS/FOM)	Has PD3 Pattern	Hes Property Sheet
Not Detabase Comments	E50:	Ves Dio	Yes No
Not Database Comments Not Contains W	E50: 8 Iarda ▼ []⊃]	Ves No	Yes No
Not Database Comments Not Contains W	ESO: s urds • []] [] Z Search]] 2 Show R	Results Dindock Page Reset Page	Yes No
Not Database Contentation Not Contents W	Sorris a second	Resulti Dudock Page Reset Page	Yes ito
Nat	BSC series > 2 2 Search > 2 How R	Next In Index Page Reset Page	Test Al.
hit	ESD 5 srds • • • 2 Search 2 Show R	Residit Dublock Page Reset Page	Type the
Hat	BSD: Sards • • •	tealt Diddd Page Read Page	The initial sector of the sect

More Composition

Chemical Formula								
Not Contains Elements •								
Empirical Formula								
Not Contains Elements 🕶		_	_		_	_		
Structural Formula								
Not Contains Elements 💌		_						
Formula Type (ANX)								
Not Contains Phrase								
Composition (1)	Numbe	r of I	Elen	nents (#	t i			
- Alas		_						
_ not	Not	2	٠	Or Or				
Weight % Atomic %	Not.	1 2 3	î	O Or				
Weight % Atomic % EL Value ESD	Not	1 2 3 4		O Or				
Weight % Atomic % EL Value ESD	Not	1 2 3 4 5 6		(a) Or				
Weight % Atomic % EL Value ESD	Not	1234567	-	O Or				

Subfiles, Quality

Search Show Results

hies/Database Hiters Periodic Table	Elements Names References Structures	Miscellaneous
atabase	Status	Ambient/Non-ambient (Amb.)
Not ICDD (00) Or ICSD-FIZ (01) Cambridge (0 NIST (02) LPF (04)	Primary Atternate Include Deleted Potterns	Net Ambient Pressure (Non-ambient) Temperature (Non-ambient) Pressure & Temperature (Non-ambi
uality Mark (QM)	Subfile/Subclass	
Net Starkeld (R) Network (D) Calculated (C) Stark (B) Calculated (C) Stark (B) Cover/Precision (D) Prototyping (P) rypothetical (r) Good (C) Minimal Acceptable	Hett Exciple Explori Foreau International Metals Micro A Micro Mic	nt ve Garage And Garag

C Undock Page Reset Page Reset All

Search V Show Results Ondock Page Reset Page Reset All

Structures Percens Symbol Code (Autorice Cargetal Data Reduced Cell Percenson Symbol Code (Percenson) The Contains Prazer (With Hydrogen) Autor-Cherded Space Croug (Nich) Autor-Cherded Space (Nich

Point and click Drop down Menus Boolean operators

6

6

IMA, IZA -Classifications and nomenclature

Commond lians		111
Not Contains Words -		
Common Name		
Not Contains Words • 🗇	6	
Mineral Name		
Not Contains Words +	6	
All Nomes		
Not Contains Words + 3		
Zeolite Classification	Mineral Classification	
Her AW - L-A(W) ▲ ○ 0: All: - APO-1 A > > ○ 0: All: - APO-1 A > > > ○ 0: All: - APO-1 A > <t< th=""><th>Her: All-Anargunstr (Grav) Could All Anargunstr (Grav) All - Andres (Grav) All - Andr</th><th></th></t<>	Her: All-Anargunstr (Grav) Could All Anargunstr (Grav) All - Andres (Grav) All - Andr	
	Structure: Please select a class	
Bilbiograp	ohic extenses	
Not Contains Words		1
Coden	Journal	
Not Exactly	Not Contains Words	
Title		
Not Contains Words		.
Chemical Abstracts Registry Service (CAS)	Year Volume	
Not Contains Phrase	Not ESD: Not Exactly -	1
Search Show Resu	ults 🎦 Undock Page Reset Page Reset All	
	14	

Structural Information in PDF-4+

Quantitative Analysis

Simulations





Basic Mineral Content in the PDF

Minerals and Related Materials – 39,410 Entries all entries have digital patterns



92 experimental PD3 patterns – mostly clays

Case Studies

Comparison of experimental data to 4 calculated data sets from four published structures of montmorillonite

* The reference publications for the structures of montmorillonite detail extensive treatments used to induce crystallinity and refine the structure. It is the opinion of many experts that these results are not accurate representations of the mineral.

Method – Analysis of a 3 phase sample of montmorillonite, from Northern China

 Match the montmorillonite
 Use the best experimental digital pattern for montmorillonite – 19 choices

Use calculated patterns, adjusted for crystallite size, for quartz and cristobalite

Scaled components - separated

– Ca0.2 (Al , Mg)2 Si4 O10 (O H)2 🗴 H2 O - 00-058-2007 (PD3, Intensity: 50.0%) — Si O2 - 04-007-2134 (Calc, Intensity: 20.0%) — Si O2 - 01-085-1780 (Calc, Intensity: 10.0%)

* A scatter function was used to adjust background levels. A more accurate representation could be made using a 2D scatter simulation for montmorillonite. X-ray diffraction from two dimensional crystals, such as montmorillonite and other minerals in the smectite group, is described by Brindley and Brown, (1980) "*Crystal Structures of Clay Minerals and their X-ray Identification*", published by The Mineralogical Society, London, UK. The program DIFFax can simulate layer stacking. http://www.public.asu.edu/~mtreacy/DIFFaX.html

Method Details

Match the best montmorillonite – use a similarity index to compare experimental data to all references

() () () () () () () () () () () () () () (11111				11111.	() () () () (111111	111111	11111	111111
Normalized R-index	PDF #	QM	Mineral Name	Chemical Formula	Compound Name	Coords	D1	D2	D3	SYS	Year
2 1.12 (2.60° - 70.00°)	00-058-2038	В	Montmorillonite, calcian	(Ca, Na) 0.3 Al2 (Si, AI) 4 O10 (Calcium Sodium Aluminum Silicate H		15.445900	4.483360	4.115850	М	1983
2 1.16 (2.60° - 70.00°)	00-058-2007	В	Montmorillonite-calcian	Ca0.2 (AI, Mg)2 Si4 O10 (OH)2 ·	Calcium Magnesium Aluminum Silicat		15.432400	4.475070	1.496970	н	1983
2.60° - 70.00°) 1.22	00-060-0318	В	Montmorillonite	(Ca, Na) 0.3 Al2 (Si, AI) 4 O10 (Sodium Calcium Aluminum Silicate H		14.024500	4.475970	2.491340	М	1983
2.60° - 56.88°) (2.60° - 56.88°)	00-060-0319	В	Montmorillonite, oriented	(Ca, Na) 0.3 Al2 (Si, AI) 4 O10 (Sodium Calcium Aluminum Silicate H		15.628900	3.103580	5.140480	M	1983
2.60° - 55.64°) (2.60° - 55.64°)	00-058-2009	В	Montmorillonite-calcian	Ca0.2 (AI, Mg)2 Si4 O10 (OH)2 ·	Calcium Magnesium Aluminum Silicat		15.673300	3.173450	5.258910	н	1983
2.60° - 57.96°) 1.38 (2.60° - 57.96°)	00-058-2008	В	Montmorillonite-calcian	Ca0.2 (AI, Mg)2 Si4 O10 (OH)2 ·	Calcium Magnesium Aluminum Silicat		15.853200	3.068780	5.149970	н	1983
2 1.44 (2.60° - 70.00°)	00-058-2010	В	Montmorillonite	Na0.3 (AI, Mg)2 Si4 O10 (OH)2 ·	Sodium Aluminum Magnesium Silicat		12.522800	3.118450	4.452390	0	1983
2.60° - 70.00°) (2.60° - 70.00°)	00-060-0315	В	Montmorillonite-Na	(Na, Ca) 0.3 Al2 (Si, AI) 4 O10 (Sodium Calcium Aluminum Silicate H		3.148310	12.478600	4.465480	M	1983
🎢 1.62 (2.60° - 58.14°)	00-060-0316	В	Montmorillonite, gly colated	(Na, Ca) 0.3 Al2 (Si, AI) 4 O10 (Sodium Calcium Aluminum Silicate H		17.108600	3.385200	8.545640	M	1983
1.62 (2.60° - 58.20°)	00-060-0320	1	Montmorillonite, gly colated, oriented	(Ca, Na) 0.3 Al2 (Si, AI) 4 O10 (Sodium Calcium Aluminum Silicate H		17.003200	3.399160	8.560500	M	1983
🎢 1.65 (2.60° - 41.92°)	00-058-2041	В	Montmorillonite, calcian, heated	(Ca, Na) 0.3 Al2 (Si, AI) 4 O10 (Calcium Sodium Aluminum Silicate H		18.521700	4.476190	7.301200	н	1983
1.74 (2.60° - 53.12°)	00-058-2011	В	Montmorillonite	Na0.3 (AI, Mg) 2 Si4 O10 (OH) 2 ·	Sodium Aluminum Magnesium Silicat		12.803800	3.164290	6.232400	0	1983
2.60° - 53.34°) (2.60° - 53.34°)	00-058-2039	В	Montmorillonite, calcian	(Ca, Na) 0.3 Al2 (Si, AI) 4 O10 (Calcium Sodium Aluminum Silicate H		12.722800	3.120050	6.246850	M	1983
2.60° - 42.82°) 1.83	00-058-2012	В	Montmorillonite, gly colated	Na0.3 (AI, Mg) 2 Si4 O10 (OH) 2 ·	Sodium Aluminum Magnesium Silicat		21.928900	4.854450	7.952340	0	1983
2 1.82 (2.60° - 42.74°)	00-058-2040	в	Montmorillonite, calcian, gly colated	(Ca, Na) 0.3 Al2 (Si, AI) 4 O10 (Calcium Sodium Aluminum Silicate H		21.633500	4.874830	7.940240	M	1983
2.60° - 42.72°)	00-058-2029	В	Montmorillonite, calcian, gly colated	Ca0.2 (AI, Mg)2 Si4 O10 (OH)2	Calcium Magnesium Aluminum Silicat		21.517500	4.888960	7.932430	х	1983
2.60° - 57.32°)	00-060-0321	в	Montmorillonite, heated, oriented	(Ca, Na) 0.3 Al2 (Si, AI) 4 O10 (Sodium Calcium Aluminum Silicate H		9.737700	3.192940	4.808500	M	1983
2.60° - 57.30°)	00-058-2013	В	Montmorillonite, heated	Na0.3 (AI, Mg) 2 Si4 O10 (OH) 2 ·	Sodium Aluminum Magnesium Silicat		9.803470	3.188360	4.843400	0	1983
2.60° - 57.36°) (2.60° - 57.36°)	00-060-0317	1	Montmorillonite, heated, oriented	(Na, Ca) 0.3 Al2 (Si, AI) 4 O10 (Sodium Calcium Aluminum Silicate H		3.193730	9.624450	4.790210	М	1983

Match

quality

Library of References

Similarity Index Matches Note: 1st choice also matched cristobalite impurity so we used the second choice without the impurity

Mine Specimens

Phase identification

Mineral Specimen from Xinjiang Gold Mine Data courtesy of Prof. Jiaju, Beijing, China

— Fe0.91 S - 04-007-0975 (Calc, Intensity: 58.3%) — (K, Na) (Fe, Al, Mg)2 (Si, Al) 4 O10 (OH)2 - 00-058-2024 (Exp-based, Intensity: 17.3%)

— Mg5 AI (Si3 AI) O10 (O H)8 - 00-029-0853 (Exp-based, Intensity: 15.0%) — Fe Ni S2 - 01-075-0611 (Calc, Intensity: 10.0%)

Plus cordierite and albite

Kimberlite Rock - Diamond bearing mineral

- Mg3 Si2 O5 (O H)4 - 00-062-0393 (Exp-based, Intensity: 55.1%) - Summation

Fosterite Magnetite (Mg) Lizardite

Kimberlite rock – 6 phases with apatite and montmorillonite

Isolated/Separated Mineral Fraction from a Kimberlite pipe core sample (clays and light density minerals removed by chemical and physical separation)

- Ca0.2 Mg1.6 Fe1.2 Al2 (Si O4)3 - 04-012-4953 (Calc, Intensity: 94.2%) — Ti O2 - 04-008-7848 (Calc, Intensity: 29.3%) — Mn3 Al2 (Si O4)3 - 04-013-2246 (Calc, Intensity: 41.2%) - Zr (Si O4) - 01-089-2655 (Calc, Intensity: 85.4%) — Fe2.939 O4 - 01-086-1356 (Calc, Intensity: 77.8%)

Kimberlite

- Indicator minerals were initially used to identify pipe locations
- Indicator minerals isolated from ore body using physical and chemical separation
- Cr analyzed in the pyrope ruby garnet
- Diamonds identified in the pipe test cores
- 20 Kimberlite pipes eventually found in the area
- None have been commercially mined

Hunt for diamonds

At present, 15 companies have diamond exploration leases for 130,000 acres of state-owned land in a 5,000-square-mile area of Iron, Dickinson and Gogebic Counties in the western section of the Upper Peninsula. An estimated 1 million acres of privately owned land also are under lease.

Chicago Tribune, June 30th 1985

TORONTO (<u>ResourceInvestor.com</u>) -- The latest results from Dianor Resources [TSX-V:<u>DOR</u>] continue to point to the company's Leadbetter property near Wawa, Ontario, Canada hosting a significantly large body of diamondiferous rock.

The company optioned an 80% stake in the road-accessible Leadbetter property in December 2004 for close to \$4 million in cash and shares over a four-year period. In excess of 330 diamonds, ranging up to 0.11 carats in size, had been collected from bedrock by prospector Joe Leadbetter in 2004. In addition, three macro diamonds were recovered from adjacent creek gravels, including a gem quality 1.39-carat stone valued at US\$3,000. A pink 0.11-carat diamond is the largest stone recovered as of April 2005 from the bedrock.

Dianor started collecting and processing samples in January of this year. Initial results caused the shares to shoot up to \$0.70 and were promising enough for the company to acquire a 70% interest in 49 claims located to the east and south of the original property, forming a contiguous land holding of approximately 16 square kilometres.

Resource Investor, May 30th, 2005 Dianor got \$30 MM in investment funding in 2010^s

Role of XRD in diamond exploration and discovery

Search for indicator minerals
 Confirmation of ruby garnets - pyrope
 Confirmation of diamonds

ADVANCED MINERAL ANALYSES BY X-RAY POWDER DIFFRACTION

- Ability to identify both crystalline and noncrystalline minerals
- Routinely identify and quantify complex multiphase specimens (2 8-phase, 2 6-phase, a 5 and a 4 phase solution in this study)
- Detection depends on resolution and counting statistics – there are both improved optics and detectors in recent years significantly lowering detection limits
- New software tools (similarity indices, whole pattern methods, multivariant analyses) greatly facilitate search, identification and quantification.