

PAUL SCHERRER INSTITUT



Excelsus
Structural Solutions

WIR SCHAFFEN WISSEN – HEUTE FÜR MORGEN



M. L. Reinle-Schmitt^a, P. P. Mazzeo^b, I. C. Madsen^c, A. Prodi^d, F. Gozzo^{b,d}

^aPaul Scherrer Institute, 5232 Villigen, Switzerland, after 1st July: Excelsus Structural Solutions (Swiss) AG, Switzerland

^bExcelsus Structural Solutions (Swiss) AG, Switzerland

^cCSIRO Process Science and Engineering, Bayview Avenue, Box 312, Clayton South 3168, Victoria, Australia

^dExcelsus Structural Solutions SPRL, Belgium

Absolute quantification of pharmaceuticals: The search of suitable internal standards

PPXRD 14, Fort Myers, Florida, 9th June 2016

This document was presented at PPXRD - Pharmaceutical Powder X-ray Diffraction Symposium

Sponsored by The International Centre for Diffraction Data

This presentation is provided by the International Centre for Diffraction Data in cooperation with the authors and presenters of the PPXRD symposia for the express purpose of educating the scientific community.

All copyrights for the presentation are retained by the original authors.

The ICDD has received permission from the authors to post this material on our website and make the material available for viewing. Usage is restricted for the purposes of education and scientific research.



PPXRD Website – www.icdd.com/ppxrd

ICDD Website - www.icdd.com

- QPA on absolute scale of traces in organic mixture
- Internal standard method

- Constraints for the choice of an internal standard:
 - Internal standard method in itself
 - Synchrotron radiation

Which strategy for absolute QPA via SR-XRPD?
Did we find **THE** internal standard ?




Photo: Paul Scherrer Institute


Outlook

1. Internal standard for pharmaceuticals and SR-XRPD
2. Search of the internal standard
3. Challenges/Preliminary results

On going project!

1. Internal standard for pharmaceuticals and SR-XRPD

- Small volume analyzed




VS




1. Internal standard for pharmaceuticals and SR-XRPD

- Small volume analyzed
- Restrictions on capillary diameter





1. Internal standard for pharmaceuticals and SR-XRPD

- Capillaries : transmission, angular resolution, modelling, highly potent and/or reactive
- Fast measurement, high angular (FWHM) resolution/tunability





- Much higher S/N and S/B
- LoQ < 0.05 wt%, LoD < 0.01 wt%



1. Internal standard for pharmaceuticals and SR-XRPD

- Poor scattering power (H, B, C, N, Si, P, S, O)
- Radiation sensitive
- Low absorption
- Low symmetry



1. Internal standard for pharmaceuticals and SR-XRPD



- Addition of a known weight ratio of internal standard to the mixture to analyze – known reference to scale all phases
- Use of intensity ratios:

SINGLE PEAK METHOD

$$w'_\alpha = \frac{I_{(hkl),\alpha} K_{(hkl),S}}{I_{(hkl),S} K_{(hkl),\alpha}} \cdot \frac{\rho_\alpha}{\rho_S} \cdot w_S$$

RIETVELT METHOD

$$w'_\alpha = \frac{S_\alpha (ZMV)_\alpha}{S_S (ZMV)_S} \cdot w_S$$



1. Internal standard for pharmaceuticals and SR-XRPD

- Direct correction for instrumental effects
- Comparable matrix effects
- Unknown compounds
- Amorphous quantification, Absolute scale
- Universal



- Internal standard tailored to analyte
- Time consuming powder processing
- Powder samples only
- Analyte mixture contamination





1. Internal standard for pharmaceuticals and SR-XRPD

• Constraints/requirements:

- High symmetry, high crystallinity
- Well-known structure
- Well known amorphous component
- Comparable dimensions

- Small isometric particle, density
- Similar mass absorption coefficient
- Small crystallite size
- Comparative intensities (wt %)

- Stable, not reactive
- Affordable, easily available



| Crystallite diameter (μm) | 40 | 10 | 1 |
|--|--------------------|--------------------|-----------------------|
| Crystallites (20 nm^3) | 5.97×10^5 | 3.82×10^7 | 3.82×10^{10} |
| Number diffracting | 12 | 760 | 38 000 |
| σ_{PS} | 0.289 | 0.036 | 0.005 |

Robert E. Dinnebier and Simon J. L. Billinge, Print ISBN: 978-0-85404-231-9, DOI:10.1039/978184755823 adapted from: D. K. Smith, Adv. X-Ray Anal., 1992, 35, 1-15

1. Internal standard for pharmaceuticals and SR-XRPD

‘Most severe limitation to QPA accuracy is particle statistics’


Deane K. Smith Powder Diffraction, 16, pp 186-191, (2001), doi:10.1154/1.1423285.

- Individual volumes representative of general phase distribution in the whole sample requires
 - Homogeneity of the analyte
 - Homogeneity of the mixture with Internal Standard
 - Consistent sampling
- Measure entire sample??



2. Search of internal standard: Analyte mixture

- QPA on the analyte mixture with internal standard candidates
- Demonstrate a method: negligible impact of test mixture
- Acceptable peak overlaps
- Particle/crystallite size
- QPA of a ternary organic mixture: APIs well-known structure
 - Majority phase: Acetaminophen (Ball milled) – 75 to 96 %w/w
 - Medium phase: Salicylic Acid (Ground) – 3 to 20 %w/w
 - Minority phase: Carbamazepine (Ground) - 0.1 to 5 %w/w



2. Search of internal standard: Literature

- **Organic?**
- **Al_2O_3**
 - Naproxen, Danazol, 10% and 50% wt mixing, Lab-machine
 - Urinary stone constituents =, 30-40 %wt mixing, Lab-machine, failure
- **LiF**
 - Sodium salt amorphous vs crystalline, 5%
 - Tolnaftate in microspheres, 20%
 - Racemic compound of ibuprofen, 20% wt mixing
- **Diamond powder, $\text{Ca}(\text{OH})_2$, $\text{CaSO}_4 \cdot 2\text{H}_2\text{O}$, ZnO, Silicon powder...**

2. Search of internal standard: 1st screening

| Internal Standard | Chemical formula | Density (g.cm ⁻³) | Crystal structure |
|---|--|-------------------------------|--|
| Hexamethylenetetramine | (CH ₂) ₆ N ₄ | 1,33 | Cubic |
| Diamond | C | 3,51 | Cubic |
| Lithium fluoride (precipitated, 99,995%, Sigma Aldrich) | LiF | 2,635 | Cubic |
| Monosodium citrate (Jungbunzlauer) | NaH ₂ C ₆ H ₅ O ₇ | 1,5 | Two known polymorphs: monoclinic & orthorombic |
| Sodium carbonate (anhydrous, ≥99,9999%, Sigma Aldrich) | Na ₂ CO ₃ | 2,54 | monoclinic or orthorombic |
| Zeolite (Faujasite) | [Na _{28.8} Ca _{14.4} (H ₂ O) ₂₆₃] [Si _{134.4} Al _{57.6} O ₃₈₄] | 1,93 | Cubic |
| Corundum (as a reference) | Al ₂ O ₃ | 4,02 | Trigonal-hexagonal |

2. Search of internal standard: 1st screening

| Internal Standard | Chemical formula | Density (g.cm ⁻³) | Crystal structure |
|---|--|-------------------------------|--------------------|
| Hexamethylenetetramine | (CH ₂) ₆ N ₄ | 1,33 | Cubic |
| Diamond | C | 3,51 | Cubic |
| Lithium fluoride (precipitated, 99,995%, Sigma Aldrich) | LiF | 2,635 | Cubic |
| Monosodium citrate (Jungbunzlauer) | | | |
| Sodium carbonate (anhydrous, ≥99,9999%, Sigma Aldrich) | | | |
| Zeolite (Faujasite) | [Na _{28.8} Ca _{14.4} (H ₂ O) ₂₆₃] [Si _{134.4} Al _{57.6} O ₃₈₄] | 1,93 | Cubic |
| Corundum (as a reference) | Al ₂ O ₃ | 4,02 | Trigonal-hexagonal |

Selected

Very hard



Selected

Low symmetry

Low symmetry

Later stage

Selected



2. Search of internal standard


- Particle size distribution: 1-5 microns
- Dry microsieving 5 to 20 microns:
 - As received, ground and ball-milled HMTA and LiF
 - agglomeration or too large *psd*
- Ball-milling
- Characterization techniques:
 - Laser granulometry
 - Optical microscope
 - Scanning electron microscopy



2. Search of internal standard

APIs

Internal standard



Amorphization ?



Electrostaticity ?
Agglomeration ?
Water uptake ?





CONSISTENT SAMPLE PREPARATION !




3. Preliminary results: HMTA

- Organic
 - similar density
 - similar mass absorption coefficient
- Cubic symmetry, highly crystalline




- Resistant to milling
 - reduced particle size distribution
 - reduced crystallite size
 - amorphization: in progress



3. Preliminary results: HMTA

- Not stable in mixture



HMTA + SA

HMTA + SA^{-Q}

LAG in EtOH



LAG in EtOH (50 μ L)

10 min grinding in agate mortar




The salt was not pure
[reactant as contaminant]

3. Preliminary results: Alumina

- NIST Quantitative Analysis Powder Diffraction Standard
- Few intense peaks






- Known structural model
- Mean particle size below 2 microns




- Well-known amorphous content (99,2% crystalline)

3. Preliminary results: Alumina

- High density, inhomogeneity along capillary






- Absorption contrast, high mass attenuation coefficient:




- Previous experiment: effect of microabsorption with 20% Al_2O_3

3. Preliminary results: LiF

- Excellent structural model (special positions)
- No preferential orientation








- Resistant to milling
 - Particle size distribution < 5 microns
 - reduced crystallite size
 - amorphization: in progress



3. Preliminary results: LiF

- Inhomogeneous distribution in spite of milling and careful mixing
- Strong agglomeration, hygroscopic







3. Preliminary results: LiF

- Confirm inhomogeneity of LiF
- Improving the statistics, insufficient volumes collected

$$Relative\ bias\ (wt\%) = \frac{weighted - refined}{weighted}$$





Relative bias wrt calculated (Rietvelt) weight ratios for LiF



3. Preliminary results: LiF







- LiF overestimated



3. Preliminary results: LiF


- PX
- ◆ ALL
- MS
- ▲ LS

- LiF overestimated
- ACP underestimated , CBZ and SA overestimated




3. Preliminary results: Analyte mixture

- Preferential orientation in Acetaminophen and Salicylic Acid



3. Preliminary results: DoC

- Degree of crystallinity method – not yet convincing results
- Problem with modelling and correction of extrinsic background



Difference > expected 0.8%

- Additional techniques:
 - Water sorption (gravimetric method)
 - Infrared spectroscopy
 - Calorimetry

Analyte mixture:



- Replace ACP, SA

Internal standard

- Identify when HMTA can be used
- Reduce agglomeration in LiF
- Amorphous

Measurements @ SLS

- Determine statistically significant number of volumes
- Compare with total powder result
- Automation of measurement
- Appropriate acquisition time
- Cluster analysis of data



- Deane K. Smith Powder Diffraction, 16, pp 186-191, (2001), doi:10.1154/1.1423285
- Wall, C. et al, Powder Technology, 09/2014; 264:409–417.
- J. W. Shell, Journal of pharmaceutical sciences 52, 1 (1963).
- Schreyer, M. et al., Journal of applied crystallography, 44, 1, 17 (2011).
- Pederson, B. M. et al., Adv. X Ray. Anal. 47, 200 (2004)
- Wandt, M.A.E. and Rodgers, A.L., *Clin. Chem.* **34/2**, 289 (1988)
- X-ray powder diffractometry, Suryanarayanan R., (1995)
- Stephenson, G. A. et al., *The Rigaku Journal*, **22**, 1 (2005)
- Dash, A.K. et al., *Journal of pharmaceutical sciences*, **91**, 4 (2002).
- Phadnis, N. V. et al., *Pharmaceutical Research*, **14**, 9 (1997)
- Clas, S. D. et al., *International Journal of Pharmaceutics*, **121**, 73 (1995)
- Rogers, T. L. et al., *Pharmaceutical Research*, **21**, 11 (2004)
- Le Troedec, M. et al., *Composites: Part A*, **39**, 514 (2008)
- Alexander L. and Klug H. P., *Anal. Chem.*, **20**, 886 (1948)
- Mandile, J. A. et al., *International Journal of Coal Geology*, **28**, 51(1995)
- Chrzanowski, F.A. et al., *Journal of pharmaceutical sciences*, **73**, 10 (1984)
- Otsuka, M. and Kaneniwa, N., *Chem. Pharm. Bull.* **31**(12), 4489 (1983)
- Riello P. et al., *J. Appl. Cryst.* 28, 121-126 (1995)

A Practical Guide for the Preparation of Specimens for X-Ray Fluorescence and X-Ray Diffraction Analysis

Victor E. Buhrke (Editor), Ron Jenkins (Editor), Deane K. Smith (Editor) , ISBN: 978-0-471-19458-3

Powder Diffraction: Theory and Practice, 2008

Robert E. Dinnebier and Simon J. L. Billinge, Print ISBN: 978-0-85404-231-9, DOI:10.1039/978184755823

Quantitative X-Ray Diffractometry, 1995

Lev S. Zevin, Giora Kimmel, ISBN: 978-1-4613-9537-9 (Print) 978-1-4613-9535-5 (Online)

X-Ray Diffraction Procedures: For Polycrystalline and Amorphous Materials, 2nd Edition, 1974

Harold P. Klug, Leroy E. Alexander, ISBN: 978-0-471-49369-3

My thanks go to

- N. Casati and MS beamline staff,
- A. Testino
- A. Weber
- E. Pomjakushina
- C. Padeste

From:
Paul Scherrer Institute





Thank you for your attention!

Contact: m.reinle-schmitt@excels.us