

WIR SCHAFFEN WISSEN – HEUTE FÜR MORGEN



Mathilde Reinle-Schmitt :: Paul Scherrer Institut

1st July 2016: Excelsus Structural Solutions (Swiss) AG

Challenges of absolute quantification of pharmaceuticals by the internal standard methods

Quantitative Phase Analysis by XRPD Workshop, PPXRD 14, Fort Myers, 6th June 2016

This document was presented at PPXRD - Pharmaceutical Powder X-ray Diffraction Symposium

Sponsored by The International Centre for Diffraction Data

This presentation is provided by the International Centre for Diffraction Data in cooperation with the authors and presenters of the PPXRD symposia for the express purpose of educating the scientific community.

All copyrights for the presentation are retained by the original authors.

The ICDD has received permission from the authors to post this material on our website and make the material available for viewing. Usage is restricted for the purposes of education and scientific research.



PPXRD Website – www.icdd.com/ppxrd

ICDD Website - www.icdd.com



Introduction

- QPA: high degree of accuracy and precision
- Focus on the internal standard method
- Appropriate choices of primary importance

Why an internal standard?
Recipe to choose it?
Case study



An internal standard for pharmaceuticals

Phase α , μ_α





An internal standard for pharmaceuticals

Phase α , μ_α



$$I_{(hkl),\alpha} = \frac{K_{1,\alpha}}{\mu_\alpha}$$

Phase β , μ_β




$$I_{(hkl),\beta} = \frac{K_{1,\beta}}{\mu_\beta}$$

An internal standard for pharmaceuticals

- Mixture

$$I_{(hkl),\alpha} = \frac{K_{1,\alpha}}{\mu_m} c_\alpha$$




$$I_{(hkl),\beta} = \frac{K_{1,\beta}}{\mu_m} c_\beta$$

An internal standard for pharmaceuticals

- Mixture


$$I_{(hkl),\alpha} = \frac{K_{1,\alpha} w_\alpha}{\rho_\alpha \mu_m^*}$$

$\mu_\alpha, c_\alpha, w_\alpha$



μ_m

$\mu_\beta, c_\beta, w_\beta$



$$I_{(hkl),\beta} = \frac{K_{1,\beta} w_\beta}{\rho_\beta \mu_m^*}$$


An internal standard for pharmaceuticals

- Mixture

$$w_\alpha = I_{(hkl),\alpha} \frac{\rho_\alpha \mu_m^*}{K_{1,\alpha}}$$

known

unknown




$$w_\beta = I_{(hkl),\beta} \frac{\rho_\beta \mu_m^*}{K_{1,\beta}}$$

An internal standard for pharmaceuticals

- Ratio


$$w_\alpha = I_{(hkl),\alpha} \frac{\rho_\alpha \mu_m^*}{K_{1,\alpha}}$$

$\mu_\alpha, c_\alpha, w_\alpha$



μ_m


$\mu_\beta, c_\beta, w_\beta$



$$w_\beta = I_{(hkl),\beta} \frac{\rho_\beta \mu_m^*}{K_{1,\beta}}$$


An internal standard for pharmaceuticals

- Ratio



An internal standard for pharmaceuticals



- Add the internal standard



$$\frac{w'_\alpha}{w_\beta} = \frac{\frac{I(hkl),\alpha}{K_{1,\alpha}} \frac{\rho_\alpha}{K_{1,\alpha}}}{\frac{I(hkl),\beta}{K_{1,\beta}} \frac{\rho_b}{K_{1,\beta}}}$$

An internal standard for pharmaceuticals



- Add the internal standard



$$\frac{w'_\alpha}{w_S} = \frac{I_{(hkl),\alpha} \frac{\rho_\alpha}{K_{1,\alpha}}}{I_{(hkl),S} \frac{\rho_S}{K_{1,S}}}$$

An internal standard for pharmaceuticals


- Add the internal standard



$$w_{\alpha}' = \frac{I_{(hkl),\alpha}}{I_{(hkl),S}} \cdot \frac{K_{1,S}}{K_{1,\alpha}} \cdot \frac{\rho_{\alpha}}{\rho_S} w_S$$

Calibration constant
from ad-hoc mixtures
with known w_S and w'_{α} :

$$\frac{K_{1,S} \rho_{\alpha}}{K_{1,\alpha} \rho_S} = \frac{I_{(hkl),S}}{I_{(hkl),\alpha}} \frac{w_S}{w'_{\alpha}}$$



An internal standard for pharmaceuticals

- Single peak method

$$w_\alpha = I_{(hkl),\alpha} \frac{\rho_\alpha \mu^*_m}{K_{1,\alpha}}$$

$$w'_\alpha = \frac{I_{(hkl),\alpha} K_{1,S}}{I_{(hkl),S} K_{1,\alpha}} \cdot \frac{\rho_\alpha}{\rho_S} \cdot w_S$$

- Whole pattern method

$$w_\alpha = S_\alpha \frac{(ZMV)_\alpha \mu^*_m}{K}$$

$$w'_\alpha = \frac{S_\alpha (ZMV)_\alpha}{S_S (ZMV)_S} \cdot w_S$$



- Absolute scale

An internal standard for pharmaceuticals

- Unknown compounds
- Rietveld: only crystalline phases

$$\sum_{i=1}^n w_i = 1$$

Without Internal Standard




WRONG!

An internal standard for pharmaceuticals

- Unknown compounds
- Amorphous quantification, Absolute scale

With Internal Standard



New weight ratios!

100 %w/w	
Amorphous, 15 w%/w	⇒ not detected
Crystalline, 22,5 %w/w	26,47 %w/w
Crystalline, 37,5 %w/w	44,11 %w/w
Crystalline, 25 %w/w	29,41 %w/w

$$\begin{aligned} &\Rightarrow 100-(22,5+37,5+25) = 15 \% \text{w/w} \\ &\Rightarrow 22,5 \% \text{w/w} \\ &\Rightarrow 37,5 \% \text{w/w} \\ &\Rightarrow 25 \% \text{w/w} \end{aligned}$$

100 %w/w

⇒ 20 w%/w
⇒ 30 %w/w
⇒ 50 %w/w

CORRECT!


C=0.85

(Correction factor: $0,25 = 0,2941 \cdot C$)

An internal standard for pharmaceuticals


- Example: absolute QPA of the following mixture:

	Phase α (minority)	Phase β	Phase γ (majority)	Internal standard
Weighted	0.09 wt%	2.38 wt%	77.527 wt%	19.997 wt%
Refined (Rietveld)	0.17 wt%	2.77wt%	76.99 wt%	20.07 wt%
Corrected	0.169 wt%	2.76wt%	76.73 wt%	19.997 wt%



An internal standard for pharmaceuticals

- Example: absolute QPA of the following mixture:



An internal standard for pharmaceuticals

- Example: absolute QPA of the following mixture:

	Phase α (minority)	Phase β	Phase γ (majority)	Internal standard
Weighted	0.09 wt%	2.38 wt%	77.527 wt%	19.997 wt%
Refined (Rietveld)	0.17 wt%	2.77wt%	76.99 wt%	20.07 wt%
Corrected	0.169 wt%	2.76wt%	76.73 wt%	19.997 wt%

Divided
by
 $\frac{20.07}{19.997} =$
 1.0036

⇒ Amorphous=100-0.169-2.76-76.73-19.997=0.34 wt%

Phase α: `corrected_weight_percent`

0.169

Phase β: `corrected_weight_percent`

2.76

Phase γ: `corrected_weight_percent`

76.73


Internal standard: `spiked_phase_measured_weight_percent` 20

Amorphous: `weight_percent_amorphous`


0.34

Amorphous problem


- Minimum amorphous content



- Not negligible



- Amorphous/crystalline ratio, DoC ?



Difference > expected 0.8%

Cline, J. P., Von Dreele, R. B., Winburn, R., Stephens, P. W. & Filliben, J. J. (2011). Acta Cryst. A67, 357-367.

Schreyer, M. et al., Journal of applied crystallography, 44, 1, 17 (2011).

An internal standard for pharmaceuticals

- Unknown compounds
- Amorphous quantification, Absolute scale
- Direct correction for instrumental effects
- Comparable matrix effects
- Universal





- Internal standard tailored to analyte
- Time consuming powder processing
- Powder samples only
- Destructive sample preparation



What about the sample?

An internal standard for pharmaceuticals


- Poor scattering power




- Longer acquisition time
- Improved signal/noise and signal/background ratios

An internal standard for pharmaceuticals

- Low absorption



- Geometry of the experiment?



An internal standard for pharmaceuticals


- 95% of organics crystallize in 5 space groups:
 – $P2_1/c$
 – $P1$
 – $P2_12_12_1$
 – $P2_1$
 – $C2/c$...
- Low symmetry

Higher symmetry




An internal standard for pharmaceuticals

- Increasing reflection density with increasing angle, where intensity drops
- Large cell size



An internal standard for pharmaceuticals

- Radiation sensitivity



- Adapt the detector!

What about the instrument?


An internal standard for pharmaceuticals with synchrotron radiation

- Synchrotron radiation, transmission, capillaries (angular resolution/counting statistics)

- High brilliance and angular resolution
- Tunability of the wavelength
- Mythen detector (also in laboratory)




- Very small sample volume, particle statistics problem
- Larger number of positions measured, automatization



Microabsorption


- Little absorption contrast required



- Underestimated weight ratio for high absorber
- Brindley correction complicated to apply
- Rather avoid than correct:
 - Similar linear absorption coefficient
 - Tuning the wavelength
 - Particle sizes

Microabsorption

- Experience with Alumina, which amount?
- Effect of particle size on QPA accuracy




- Upper limit to the particle size

	Density (g.cm ⁻³)	Attenuation length (CXRO, in microns)
Paracetamol	1.26	5120
Salicylic Acid	1.44	4000
Carbamazepine	1.29	5880
Corundum	4.02	299

Adapted from: Pederson et al., Advances in x-ray analysis, 46, 2003.
 Pederson et al., Advances in x-ray analysis, 47, 2004.

Microstructure: Particle size



- Precision of QPA: homogeneity / particle statistics
- Closest particle size, similar density and particle shape
- No extinction effects
- Preferential orientation: isometric particles



- Preparation step mandatory
 - Grinding
 - Ball milling...

Microstructure: Crystallite size

- Avoid fluctuations along Debye-Scherrer rings



- Avoid fluctuations between replicate samples

Crystallite diameter (μm)	40	10	1
Crystallites (20 mm^3)	5.97×10^5	3.82×10^7	3.82×10^{10}
Number diffracting	12	760	38 000
σ_{PS}	0.289	0.036	0.005

$$\sigma_{PS} = \sqrt{N_{\text{diff}}}/N_{\text{diff}}$$


Symmetry

- Limit peak overlap
- High symmetry
- High crystallinity




Symmetry

- Limit peak overlap
- High symmetry
- High crystallinity




Symmetry

- Limit peak overlap
- High symmetry
- High crystallinity



- Angular range
- Well-known structure
- Known crystallinity



Additional constraints

- Stability
 - during preliminary powder preparation/mixing
 - in mixture
 - to x-rays
 - in time during storage
- Safe to handle
- Cheap, easily available
- Not present in the analyte mixture !


Formulated drugs

besides internal standard addition
analyte mixture should stay **identical!**


Amount of internal standard

- Comparable diffracted intensities

ALWAYS TEST FIRST !



Intensity too high!



More appropriate!

Case study

- QPA on absolute scale of traces in organic mixture with synchrotron radiation
- Constraints due to:
 - ✓ Internal standard method
 - ✓ Use of SR-XRPD
 - ✓ Application on pharmaceuticals

Which strategy for SR-XRPD?

Did we find **THE** internal standard ?





Photo: Paul Scherrer Institute

Case study: Analyte mixture

- QPA of a ternary organic mixture with several candidates
- Negligible impact of intrinsic properties on the refinement
- APIs well-known structure:
 - Majority phase: Acetaminophen (Ball milled) – 75 to 96 %w/w
 - Medium phase: Salicylic Acid (Ground) – 3 to 20 %w/w
 - Minority phase: Carbamazepine (Ground) - 0.1 to 5 %w/w





- Adapt particle and crystallite size

Case study: Internal standard candidates

Internal Standard	Chemical formula	Density (g.cm ⁻³)	Crystal structure
Hexamethylene-tetramine	(CH ₂) ₆ N ₄	1,33	Cubic
Diamond	C	3,51	Cubic
Lithium fluoride (precipitated, 99,995%, Sigma Aldrich)	LiF	2,635	Cubic
Monosodium citrate (Jungbunzlauer)	NaH ₂ C ₆ H ₅ O ₇	1,5	Two known polymorphs: monoclinic & orthorombic
Sodium carbonate (anhydrous, ≥99,9999%, Sigma Aldrich)	Na ₂ CO ₃	2,54	monoclinic or orthorombic
Zeolite (Faujasite)	[Na _{28.8} Ca _{14.4} (H ₂ O) ₂₆₃] [Si _{134.4} Al _{57.6} O ₃₈₄]	1,93	Cubic
Corundum (as a reference)	Al ₂ O ₃	4,02	Trigonal-hexagonal

Case study: Internal standard candidates

Internal Standard	Chemical formula	Density (g.cm ⁻³)	Crystal structure
Hexamethylene-tetramine	(CH ₂) ₆ N ₄	1,33	Cubic
Diamond	C	3,51	Cubic
Lithium fluoride (precipitated, 99,995%, Sigma Aldrich)	LiF	2,635	Cubic
Monosodium citrate (Jungbunzlauer)			
Sodium carbonate (anhydrous, ≥99,9999%, Sigma Aldrich)			
Zeolite (Faujasite)	[Na _{28.8} Ca _{14.4} (H ₂ O) ₂₆₃] [Si _{134.4} Al _{57.6} O ₃₈₄]	1,93	Cubic
Corundum (as a reference)	Al ₂ O ₃	4,02	Trigonal-hexagonal

Selected

Very hard

Selected

Low symmetry



Low symmetry

Later stage

Selected

Case study: Sample preparation


- Target particle size distribution: 1-5 microns
- Microsieving 5 to 20 microns:
 - As received, ground and ball-milled HMTA and LiF
 - agglomeration problem / too large psd
- Ball-milling
- Characterization techniques:
 - Laser granulometry
 - Optical microscopy
 - Scanning electron microscopy



Case study: Sample preparation

APIs

Internal standard



Amorphization ?





Electrostaticity ?
Agglomeration ?
Water uptake ?



**CONSISTENT SAMPLE
PREPARATION !**

Case study: SR-XRPD


- Avoid transparency effect: transmission geometry
- Achieve sufficient statistics in small volumes:
 - Trade off between capillary diameter and achieved resolution
 - Capillary spinning
 - Data on multiple cap. volumes



Case study: HMTA

Crystal structure	PO	Reactivity	Agglomerates	Particle/crystallite size	Homogeneity	Micro-absorption	QPA analysis
✓	no	⚠	no	no	no	no	?




- HMTA: not stable in mixture !



Case study: Al_2O_3

Crystal structure	PO	Reactivity	Agglomerates	Particle/crystallite size	Homogeneity	Micro-absorption	QPA analysis
✓	no	✓	no	✓	⚠	no	Al_2O_3 systematically overestimated

- Inhomogeneity along capillary





- High density

Case study: LiF

Crystal structure	PO	Reactivity	Agglomerates	Particle/crystallite size	Homogeneity	Micro-absorption	QPA analysis
✓	no	✓	⚠	✓	⚠	no	LiF systematically overestimated



- Inhomogeneous distribution in spite of milling and careful mixing



Case study: LiF

Crystal structure	PO	Reactivity	Agglomerates	Particle/crystallite size	Homogeneity	Micro-absorption	QPA analysis
✓	no	✓	⚠	✓	⚠	no	LiF systematically overestimated

- Inhomogeneous distribution in spite of milling and careful mixing
- Strong agglomeration, hygroscopic



Case study: Preliminary results

	Alternative?		
	Acetaminophen	Salicylic Acid	Carbamazepine
Crystal structure	✗	✗	✓
PO	yes	yes	no
Reactivity	✓	✓	✓
Agglomerates	no	no	no
Particle/ Crystallite size	✗	✗	✓
Homogeneous distribution	✗	✗	no
Micro-absorption	no	no	no
QPA satisfying	always underestimated	always overestimated	always overestimated

Conclusion

Constraints must be taken into account due to

- Internal standard's role in itself
- Analyte mixture
- Instrument used
 - ⇒ Adapt internal standard to analyte and instrument

Step-by-step approach, preliminary tests on important points

Know your pure phases the best you can (amorphous content)

Consistency in sample preparation

References

- Deane K. Smith Powder Diffraction, 16, pp 186-191, (2001), doi:10.1154/1.1423285
Wall, C. et al, Powder Technology, 09/2014; 264:409–417.
J. W. Shell, Journal of pharmaceutical sciences 52, 1 (1963).
Schreyer, M. et al., Journal of applied crystallography, 44, 1, 17 (2011).
Pederson, B. M. et al., Adv. X Ray. Anal. 47, 200 (2004)
Wandt, M.A.E. and Rodgers, A.L., *Clin. Chem.* **34/2**, 289 (1988)
X-ray powder diffractometry, Suryanarayanan R., 1995
Stephenson, G. A. et al., *The Rigaku Journal*, **22**, 1 (2005)
Dash, A.K. et al., *Journal of pharmaceutical sciences*, **91**, 4 (2002).
Phadnis, N. V. et al., *Pharmaceutical Research*, **14**, 9 (1997)
Clas, S. D. et al., *International Journal of Pharmaceutics*, **121**, 73 (1995)
Rogers, T. L. et al., *Pharmaceutical Research*, **21**, 11 (2004)
Le Troedec, M. et al., *Composites: Part A*, **39**, 514 (2008)
Alexander L. and Klug H. P., *Anal. Chem.*, **20**, 886 (1948)
Mandile, J. A. et al., *International Journal of Coal Geology*, **28**, 51(1995)
Chrzanowski, F.A. et al., *Journal of pharmaceutical sciences*, **73**, 10 (1984)
Otsuka, M. and Kaneniwa, N., *Chem. Pharm. Bull.* **31**(12), 4489 (1983)
Riello P. et al., *J. Appl. Cryst.* **28**, 121-126 (1995)

A Practical Guide for the Preparation of Specimens for X-Ray Fluorescence and X-Ray Diffraction Analysis
Victor E. Bahrke (Editor), Ron Jenkins (Editor), Deane K. Smith (Editor) , ISBN: 978-0-471-19458-3

Powder Diffraction: Theory and Practice, 2008
Robert E. Dinnebier and Simon J. L. Billinge, Print ISBN: 978-0-85404-231-9, DOI:10.1039/978184755823

Quantitative X-Ray Diffractometry, 1995
Lev S. Zevin, Giora Kimmel, ISBN: 978-1-4613-9537-9 (Print) 978-1-4613-9535-5 (Online)

X-Ray Diffraction Procedures: For Polycrystalline and Amorphous Materials, 2nd Edition, 1974
Harold P. Klug, Leroy E. Alexander, ISBN: 978-0-471-49369-3

My thanks go to

- Ian Madsen, CSIRO Process Science and Engineering , Australia
- N. Casati and MS beamline staff, PSI, Switzerland
- A. Testino PSI, Switzerland
- A. Weber PSI, Switzerland
- E. Pomjakushina, PSI, Switzerland
- C. Padeste, PSI, Switzerland
- P. Mazzeo, A. Prodi, F. Gozzo, Excelsus Structural Solutions



Wir schaffen Wissen – heute für morgen



Thank you for your attention!

Contact: m.reinle-schmitt@excels.us