


Sample Preparation and Laboratory Instrumentation



Arnt Kern, Pamela Whitfield



This document was presented at PPXRD - Pharmaceutical Powder X-ray Diffraction Symposium

Sponsored by The International Centre for Diffraction Data

This presentation is provided by the International Centre for Diffraction Data in cooperation with the authors and presenters of the PPXRD symposia for the express purpose of educating the scientific community.

All copyrights for the presentation are retained by the original authors.

The ICDD has received permission from the authors to post this material on our website and make the material available for viewing. Usage is restricted for the purposes of education and scientific research.



PPXRD Website – www.icdd.com/ppxrd

ICDD Website - www.icdd.com

Introduction




- Consider sampling issues
 - Sample-related problems
 - Appropriate sample preparation
 - Adequate particle statistics
 - Random particle orientation
- Think about the sample presentation (instrument) geometry
 - Bragg-Brentano geometry
 - Always focusing
 - Always flat plate reflection

Highest intensity (LLOD, LLOQ)
Best resolution
 - Debye-Scherrer geometry
 - Parallel or focusing beam
 - Flat plate reflection (parallel beam)
 - Foil or capillary transmission (parallel and focusing beam)

Minimal sample issues

Bragg-Brentano Geometry



Bragg-Brentano Geometry




- Focusing geometry, flat plate reflection
- Highest intensity, best resolution
 - Best LLODs and LLOQs
- Requires relatively large sample amounts
- Preferred orientation issues
- Sample transparency issues (organic materials!)



Top-loading






- Simplest but most prone to inducing preferential orientation
- Special holders often in this category



Alternative holders such as cavity zero background silicon or air-sensitive often top-loading as well



Back-loading




Flat Plate Reflection




For accurate peak intensities

- Use an infinitely thick sample
- High angle intensities reduced if required thickness inadequate
- Beam footprint must not exceed sample length at low angles



Thick sample: constant diffraction volume $V_2 = V_1$




Thin sample: lower diffraction intensity at higher angles $V_2 < V_1$

Flat Plate Reflection




- Fixed divergence slits




- Constant diffraction volume!

- Variable divergence slits



- Variable diffraction volume requires SW correction!
- Avoid if possible!

Debye-Scherrer Geometry



Johansson
monochromator

LYNXEYE

Motorized slit

Capillary
sample
stage

Debye-Scherrer Geometry



- Parallel beam - Göbel mirror
 - Modest resolution, good intensity
 - Flat plate reflection, foil or capillary transmission
- Focusing beam - focusing Göbel mirror
 - Intermediate resolution, high intensity
 - Foil or capillary transmission
- Focusing beam - primary monochromator
 - High resolution, poor intensity, pure $K\alpha_1$
 - Foil or capillary transmission
- Requires small sample amounts but tedious sample preparation
- Allows to minimize preferred orientation issues
- XRPD + SAXS (Cu radiation), XRPD and PDF (Mo radiation)

Foil Transmission




- Sample is sprinkled between 2 foils, and X-ray beam shoots through
- Some instruments can mix flat plate reflection and foil transmission sample holders in sample changers
- Ideal for organic and other low absorbing samples



Capillary Transmission



- XRPD + SAXS (Cu radiation), XRPD and PDF (Mo radiation)
- Most effective at reducing preferred orientation effects (except needle-like morphology)
- Made from glass, borosilicate, quartz, Kapton or other low absorbing materials
- Different sizes (0.1 – 1mm diameter) depending on absorption





Sources of Error


Preferred Orientation



- Bragg-Brentano geometry
- Parallel or focusing beam Debye-Scherrer geometry using capillaries



Rotation parallel to the scattering vector does not minimize preferred orientation effects!



www.bruker.com